WESTERN DETECTION

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










주요 이미징 어플리케이션				플레이트 이미징 어플리케이션	
					
Fluorescent western blots	Chemiluminescent western blots	Colorimetrically stained western blots	Combined fluorescent and chemiluminescent western blots	Visible colonies (e.g., stained with crystal violet dye)	Fluorescent colonies (e.g., GFP expression)
				정성 가시화 (Qualitative visible)* 이미징 어플리케이션	
					
Fluorescently stained nucleic acid gels	Fluorescently stained protein gels	Colorimetrically stained protein gels	Colorimetrically stained membranes	Leaf sections	Thin layer chromatography (TLC) plates

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[White paper: Colony counting with iBright Imaging Systems](#)

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[White paper: Optimized gel imaging and analysis with iBright Tray Adapters for E-Gel Agarose Gels](#)

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
그림 12. 형광 멀티 플렉싱으로 다른 형광 채널을 이용해 각 단백질 신호를 캡처해 stripping이나 reprobing 없이 동일 블랏에서 2개 이상의 단백질을 검출할 수 있습니다. Human colon cancer cells (HCT116) were serum-starved for 24 hours and pretreated with LY924002 (50 μ M, 1 hr), rapamycin (10 nM, 1 hr) and/or BEZ235 (500 nM, 1 hr). Following pretreatment, insulin-like growth factor-1 (hIGF-1) was added to each sample (12.8 nM, 15 min). Cells were lysed and prepared for reducing SDS-PAGE and 20 μ g of each sample was electrophoresed on a Novex 4-20% Tris-Glycine Gel, WedgeWell format. The protein was transferred to a PVDF membrane and the resulting blot was blocked with Thermo Scientific™ Blocker™ FL buffer and then probed with the following primary antibodies overnight: anti-pAKT rabbit mAb and anti- β -actin mouse mAb. The blot was washed and probed with the following secondary antibodies for one hour: Invitrogen™ Goat Anti-Rabbit IgG-Alexa Fluor™ Plus 546 and Goat Anti-Mouse IgG-Alexa Fluor™ Plus 800 conjugates. The blot was then washed and imaged on the iBright FL1500 imager using appropriate settings. The composite image shows the overlaid channels from each protein (p-AKT: Alexa Fluor Plus 800 conjugate and β -actin: Alexa Fluor Plus 546 conjugate).

그림 13. 멀티플렉싱의 활용으로 복잡한 생물학적 신호전달 경로에 대한 연구를 웨스턴 블랏으로 수행할 수 있습니다. Western blot analysis demonstrates the phosphorylation of AKT and GSK-3 $\alpha\beta$ after IGF-1 treatment of human colon cancer cells (HCT116). HCT116 cells were serum-starved for 24 hours and pretreated with the PI3K pathway inhibitor LY924002 (50 μ M, 1 hr). Following pretreatment, insulin-like growth factor-1 (hIGF-1) was added to each sample (12.8 nM, 15 min) to activate the PI3K pathway. Cells were lysed and prepared for reducing SDS-PAGE and 30 μ g of each sample was electrophoresed on a Novex 4-20% Tris-Glycine Gel, WedgeWell format. The protein was transferred to a PVDF membrane and the resulting blot was blocked with Blocker FL buffer. Blot 1 was probed overnight with primary antibodies anti-pAKT rabbit mAb and anti-AKT mouse mAb. Blot 1 was washed and probed with Invitrogen™ Goat Anti-Rabbit Alexa Fluor™ Plus 800 and Goat Anti-Mouse Alexa Fluor™ Plus 647 conjugates. Blot 1 was washed again and probed with Invitrogen™ anti- β -actin mAb DyLight™ 488 conjugate for one hour. Blot 1 was washed and imaged on the iBright FL1500 imager using appropriate settings. Blot 2 was probed overnight with primary antibodies anti-p-GSK-3 $\alpha\beta$ rabbit mAb and anti-GSK-3 $\alpha\beta$ mouse mAb. Blot 2 was washed and probed with Invitrogen™ Goat Anti-Rabbit IgG-Alexa Fluor™ Plus 647 and Goat Anti-Mouse IgG-Alexa Fluor Plus 800 conjugates. Blot 2 was washed and imaged on the iBright FL1500 imager using appropriate settings.




그림 14. 형광 멀티플렉싱은 동일 블랏에서 비슷한 분자량을 가진 여러 타겟을 깔끔하게 구별할 수 있습니다. A composite image is shown along with images showing the single-color signals of individual proteins. Visualizing the individual signals can enable assessment of details that may be harder to see in a composite. PDGF was used to induce phosphorylation of AKT (Ser 473) in mouse embryo fibroblast cells (NIH-3T3). The cells were lysed and prepared for reducing SDS-PAGE and electrophoresed on a Novex 4-20% Tris-Glycine Gel, WedgeWell format. The proteins were transferred to a PVDF membrane and the resulting blot was blocked with Blocker FL buffer and probed with the following primary antibodies overnight: anti-pAKT rabbit mAb and anti-AKT mouse mAb. The blot was washed and probed with the following secondary antibodies for one hour: Goat Anti-Rabbit IgG-Alexa Fluor Plus 647 and Goat Anti-Mouse IgG-Alexa Fluor Plus 800 conjugates. The blot was washed and imaged on the iBright FL1500 imager using appropriate settings.

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
Four-color fluorescent blot


Alexa Fluor Plus 488 dye—RB1 protein with GST and HA tag (134 kDa)

Alexa Fluor 546 dye—calreticulin protein (55 kDa)

Alexa Fluor Plus 680 dye—p23 protein (23 kDa)

Alexa Fluor Plus 800 dye—HSP90 protein (90 kDa)

Membrane—iBright prestained marker

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Excitation channel	Filter range (nm)	Emission channel	Filter range (nm)	Compatible fluorophores
EX1	455–485	EM1	508–557	Alexa Fluor 488, Alexa Fluor Plus 488 dyes
EX2	515–545	EM2	568–617	Alexa Fluor 546, Alexa Fluor Plus 555 dyes
EX3	608–632	EM3	675–720	Alexa Fluor 647, Alexa Fluor Plus 647 dyes
EX4	610–660	EM4	710–730	Alexa Fluor 680, Alexa Fluor Plus 680 dyes
EX5	745–765	EM5	800–850	Alexa Fluor 790, Alexa Fluor Plus 800 dyes
EX (Green Trans)	490–520	EM2	568–617	Ethidium bromide, SYBR Safe, SYBR Gold, SYPRO Ruby, SYPRO Orange dyes

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thermofisher.com/digitalsmartstartCL750를 방문하세요.

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thermofisher.com/abassurance방문하세요.

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	iBright CL750 Imaging System	iBright CL1500 Imaging System	iBright FL1500 Imaging System
	Essential western blot and gel imaging—transition from darkroom and film-based detection with ease	Expanded application support with many of the same high-performance specifications as the iBright FL1500 Imaging System	Maximum application support, including fluorescent western blot imaging with up to 4 fluorophores at a time
Camera			
Detector	Cooled 16-bit CCD; 65,535 shades of gray		
Resolution	9.1 megapixels		
Lens	Fixed, 25 mm, f/0.95		
Field of view	22.5 x 18.0 cm (W x D); image up to 4 mini blots or gels		
Binning modes	1 x 1, 2 x 2, 4 x 4 (high resolution and sensitivity settings)	1 x 1, 2 x 2, 3 x 3, 4 x 4, 5 x 5, 6 x 6, 8 x 8 (maximum flexibility for adjusting resolution and sensitivity)	
Zoom	1–2x, digital (digital zoom reduces the effective resolution of the zoomed image)	1–8x (1–2x mechanical, 1–4x digital); mechanical zooming improves sensitivity by moving the camera closer to the sample stage and reducing the focal length	
System interface			
Touchscreen	12.1-inch capacitive LCD display; 1,024 x 768 pixels		
Storage and connectivity			
USB	2 x USB 2.0		
Networking	Ethernet port, cloud-based connectivity; optional Wi-Fi adapter sold separately		
Image file formats	G2i (proprietary), TIFF, JPG, PNG		
Hard drive	64 GB SSD	256 GB SSD	
System hardware			
Sample drawer and stage	Manually operated drawer with fixed stage	Automatic drawer with automatic rotating sample stage	
Filter sets	2 filters (0 excitation, 2 emission)	4 filters (2 excitation, 2 emission)	12 filters (6 excitation, 6 emission)
Illumination sources	• Green LED transilluminator	• Green LED transilluminator • Epi white LED	• Green LED transilluminator • Epi white LED • Epi near-IR LED
System software			
Automated features and algorithms	• Automatic zoom • Automatic focus • Automatic exposure (Smart Exposure) • Automatic onboard image analysis • High dynamic range (HDR) image capture capability (Smart Range HDR*)		
Standalone analysis applications	iBright Analysis Software (desktop version for macOS or Windows operating system), Thermo Fisher Connect Platform (cloud-based), and iBright Analysis Software—Secure (desktop version supporting 21 CFR Part 11 compliance)		
Core imaging applications			
Colorimetric protein gel imaging	•	•	•
Fluorescent protein gel imaging	•	•	•
Fluorescent nucleic acid gel imaging	•	•	•
Colorimetric membrane imaging	Limited**	•	•
Chemiluminescent western blot imaging	•	•	•
Colorimetric western blot imaging	–	•	•
Fluorescent western blot imaging	–	–	•
Specialty plate–based imaging applications			
Fluorescent colony counting	•	•	•
Visible colony counting	•	•	•
Qualitative visible imaging applications†			
Opaque objects	–	•	•
Shipping, warranty, and upgrade			
Upgrade option	No	Upgrade to the iBright FL1500 system‡	Not applicable
Weight	Approximately 47 kg (105 lb)	Approximately 50 kg (110 lb)	
Warranty	1 year from date of purchase	2 years from date of purchase	
Dimensions (L x W x H)	68 x 38 x 60 cm		

* Smart Range HDR available in chemiluminescent blot mode only.

** Membranes must be imaged when wet using Thermo Scientific™ Pierce™ Reversible Protein Stain Kit for Nitrocellulose Membranes or Thermo Scientific™ Pierce™ Reversible Protein Stain Kit for PVDF Membranes.

† Applications enable qualitative visualization of objects or confirmation of signal. Not recommended for quantitation.

‡ Upgrades not available in all regions. Please check with your sales representative for upgrade details.

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Ordering information

Product	Cat. No.
iBright FL1500 Imaging System	
1 instrument, including SmartStart Orientation and 2-year warranty	A44115
iBright CL1500 Imaging System	
1 instrument, including SmartStart Orientation and 2-year warranty	A44114
iBright CL750 Imaging System	
1 instrument, including 1-year warranty and digital SmartStart Orientation	A44116
iBright SAE Software for 21 CFR Part 11	
1 license (one license required per instrument)	A49208
iBright Tray Adapters for E-Gel Precast Agarose Gels	
iBright Tray Adapter for E-Gel Agarose Gels, 11/22 well	A56600
iBright Tray Adapter for E-Gel Agarose Gels, 48/96 well	A56599

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