CHAYON

VANTastar

2024. 06. 27.

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End time of kinetic window 1 (s): 0.5 → $180(\text{회 반복}) \times 0.5(\text{초 간격}) / \text{well} = 90 \text{ 초/well}$

Min. interval time:

Pause before - for 0 seconds

Check timing ☒ Use enhanced dynamic range

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Fluorescence Intensity - Spectral Scan

Basic Parameters | Layout | Concentrations & Volumes | Shaking

Protocol name: 1. NEW TEST

Microplate: GREINER 384 SMALL VOLUME

2. Optic
☒ Top optic ☐ Bottom optic

3. Optic Settings
Monochromator settings:

General Settings
Settling time (0.0...1.0 s): 0.1

4. Measurement start time (0...1200.0 s): 0.0
No. of flashes (0...200): 20

Comment


☐ Pause before plate reading for 0 seconds

Check timing Start measurement OK Cancel Help

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Sample Blank Standard
Control Pos.Ctrl. Neg.Ctrl.
Empty

Groups
☐ On

Index
Start value: 97
☐ Constant ☒ Increase

Replicates
Number: 1
☒ Horizontal ☐ Vertical

Reading direction:

96	1	2	3	4	5	6	7	8	9	10	11	12
A	X1	X2	X3	X4	X5	X6	X7	X8	X9	X10	X11	X12
B	X13	X14	X15	X16	X17	X18	X19	X20	X21	X22	X23	X24
C	X25	X26	X27	X28	X29	X30	X31	X32	X33	X34	X35	X36
D	X37	X38	X39	X40	X41	X42	X43	X44	X45	X46	X47	X48
E	X49	X50	X51	X52	X53	X54	X55	X56	X57	X58	X59	X60
F	X61	X62	X63	X64	X65	X66	X67	X68	X69	X70	X71	X72
G	X73	X74	X75	X76	X77	X78	X79	X80	X81	X82	X83	X84
H	X85	X86	X87	X88	X89	X90	X91	X92	X93	X94	X95	X96

2.

Check timing ☒ Use enhanced dynamic range Start measurement **OK** Cancel Help

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Fluorescence Intensity - Plate Mode

Basic Parameters | Layout | Concentrations & Volumes | **Shaking** | Output

Shaking Actions

1.

Shaking 1

2.

Shake: Shake between readings ▾

Shake mode: Double orbital ▾

Frequency: 700 rpm ▾

Cycle(s): all cycles

Time: remaining time between cycles

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Check timing Start measurement **OK** Cancel Help

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Dynamic range: Default focal height (11 mm) New focal height...

Start Measurement - Nano DRL F

Focus and Dynamic Range / Plate IDs Sample IDs Dilution Factors Crosstalk Determination

Change layout

96	1	2	3	4	5	6	7	8	9	10	11	12
A	X1	X2	X3	X4	X5	X6	X7	X8	X9	X10	X11	X12
B	X13	X14	X15	X16	X17	X18	X19	X20	X21	X22	X23	X24
C	X25	X26	X27	X28	X29	X30	X31	X32	X33	X34	X35	X36
D	X37	X38	X39	X40	X41	X42	X43	X44	X45	X46	X47	X48
E	X49	X50	X51	X52	X53	X54	X55	X56	X57	X58	X59	X60
F	X61	X62	X63	X64	X65	X66	X67	X68	X69	X70	X71	X72
G	X73	X74	X75	X76	X77	X78	X79	X80	X81	X82	X83	X84
H	X85	X86	X87	X88	X89	X90	X91	X92	X93	X94	X95	X96

2. X59

Focus: New focal height...

Dynamic range: Enhanced dynamic range (EDR)

Monochromator / Filter Settings Gain

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Plate Identification

ID1: <protocol> ID2: <date>, <time> ID3: <1+####>

Automatically enter the plate IDs previously used with this protocol

No. of executed runs since program start: 0 Total no. of executed runs: 0

Aperture 96/384 recommended

Delay: 0 s Start measurement Save & Close Cancel Help

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
- Kinetics mode: 샘플과 같은 농도의 형광물질이 포함된 blank well 선택 후, target value를 10%로 set 합니다.
- As all the samples except blank have the same concentration of fluorophore in them at the beginning you should set gain at 10% required value and then do a gain adjustment on any well that has the fluorophore. This will set the initial value at roughly 26000 (10%) of the full range for the first reading
- As the fluorescence increases then the window will show values until it gets to 260,000 so you have a big window to work with.

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LB2 - CLARIOstar MARS

File Home View Calculations Templates Layout Formats and Settings

Open Close Define Print Report Quick Print Export Excel Report Close All Sign Change User Wizard Calculations Ranges Export Table to Excel Export Table to ASCII Manage Wavelengths Layout Add Bar Chart Templates Add Button

Navigation / Data Selection

Available Data

- LB2 (66)
 - Test Settings
 - Layout
 - Data
 - Raw Data
 - F: 450-80 [1]
 - F: LP 610 [2]
 - Calculations
 - Ratio
- LB2 (70)

Test Name: LB2
ID1: LB2; ID2: 2019-10-24, 오전 11:48:16
Luminescence

Display Mode: ☒ Values ☐ Legend in I

Excel Report Settings

Select Page for Export

- ☒ Microplate View
- ☐ Table view
- ☐ Transpose table (swap rows and columns)
- ☐ Fit results
- ☒ Protocol information
- ☐ Microplate as bitmap
- ☐ Bar chart as bitmap
- ☐ Signal curve as bitmap
- ☐ Scan curve as bitmap
- ☐ Standard fit curve as bitmap

Common Excel Export Settings

- ☐ Use one workbook for each page
- ☒ Add test run information Position: above data

Microplate View Excel Export Settings

- ☐ Export all cycles/intervals
- ☐ Export data for each cycle/interval in a separate excel sheet
- ☒ Export data for all cycles/intervals in one excel sheet (one below the other)
- ☐ Export detailed well scan data

Export report to Excel Cancel OK Help

	7	8
36	1.93	
33	5670400	75
30	10916	
31	1.52	
26	2471646	23
27	2746	

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LB2 - CLARIOstar MARS

File Home View Calculations Templates Layout Formats and Settings

Open Close Define Print Report Quick Print Export Excel Report Sign Change User Wizard Calculations Ranges Export Table to Excel Export Table to ASCII Manage Wavelengths Layout Add Bar Chart Templates Add Button

Navigation / Data Selection

Available Data

- LB2 (66)
 - Test Settings
 - Layout
 - Data
 - F: 450-80 [1]
 - F: LP 610 [2]
 - Calculations
 - Ratio
- LB2 (70)

Test Name: LB2
ID1: LB2; ID2: 2019-10-24, 오전 11:48:16
Luminescence

Define Print Report

Printer Settings

Print orientation: ☐ Landscape ☒ Portrait

Print destination: PDF File

Available Print Objects:

- Microplate View
- Table View
- Protocol Information
- 21 CFR part 11
- Comment
- Detailed Information of:
 - Blank corrected based on Raw Data (485/535)
 - Spacing
 - Horizontal line
 - Image

Selected:

- OneNote
- 복사기
- OneNote 2013으로 보내기
- Microsoft XPS Document Writer
- Microsoft Print to PDF
- HP Color LaserJet CP3525 PCL6 Class Driver
- Hewlett-Packard HP LaserJet P2035
- Hancom PDF
- Fax
- PDF File
- HTML File(s)

Page No.: 1 of 2

Preview Mode: High Quality Speed




☒ Open PDF file after creating.

Header / Footer Print Settings... Page Setup... Print pages: Save as PDF Close Help

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The screenshot shows the CHAYON software interface with the 'Manage Protocols' dialog box open. The 'Measurement Method and Mode' dialog box is also open, showing the 'Absorbance' method and 'Plate mode (slow kinetics)' reading mode selected. The 'Test Protocols' table lists various protocols, with 'Growth OD600' highlighted.

Measurement Method and Mode

Measurement Method

- ☐ Fluorescence Intensity
- ☐ Time Resolved Fluorescence
- ☐ Fluorescence Polarization
- ☐ Luminescence
- ☒ Absorbance
- ☐ AlphaScreen

Reading Mode

- ☐ Endpoint
- ☒ Plate mode (slow kinetics)
- ☐ Well mode (flash kinetics)
- ☐ Spectral scan

Test Protocols

Protocol Name	Method	Mode	Optic	Microplate
BCA	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
beta-gal	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
Bradford	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
Elisa 405	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
Elisa 450	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
Elisa 492	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
➤ Growth OD600	Absorbance	Plate mode	Spectr.	SBS STANDARD 96
Lowry	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
MTT toxicity	Absorbance	Endpoint	Spectr.	SBS STANDARD 96
NADH	Absorbance	Plate mode	Spectr.	SBS STANDARD 96
dsDNA	Absorbance Spectra	Endpoint	Spectr.	SBS STANDARD 96
LVis Holmium-Filter	Absorbance Spectra	Endpoint	Spectr.	BMG LABTECH 96
LVis ND-Filters	Absorbance Spectra	Endpoint	Spectr.	BMG LABTECH 96
Protein	Absorbance Spectra	Endpoint	Spectr.	SBS STANDARD 96

Buttons: New, Edit, Copy, Export, Import, Delete, Close, Help

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Start Measurement - Nano DRL F

Focus and Dynamic Range / Plate IDs | Sample IDs / Dilution Factors | Crosstalk Determination

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Plate Identification

ID1: <protocol> ID2: <date>, <time> ID3: <1+####>

☒ Automatically enter the plate IDs previously used with this protocol

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
RNA LVIS Test - SPECTROstar Nano MARS

Test Name: RNA LVIS Test (based on Template: RNA)

Absorbance spectrum

Wavelength: 220 nm <Select a Wavelength group>

Drag a column header here to group by that column

Manage Wavelength(s)

Automatic wavelength detection

Minimum and maximum calculations based on:

Raw Data Well: A10

from 220 nm to 360 nm

☐ at next local minimum ☐ add Wavelengths for all local minima

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Open Close Define Print Report Quick Print Export Excel Sign Wizard Calculations Ranges Export Table Export Table Manage Layout Abs in Add Bar

Navigation / Data Selection Available Data

Manage Test Runs

Recent Test Runs Microplates LVis

LVis Plate measurements (1):

Drag a column header here to group by that column

Test ID	Test Name	ID 1	ID 2	ID 3	Date	Time	Measurement Method	Signed	Protocol Comment	State
16	Protein Quant	BSA Std Curve 2			2012-01-04	오후 1:2...	Absorbance			layout ...

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Copy Delete Export Import Merge Test Runs Open