



User Manual

Milli-Q® Direct 8/16 System



About this User Manual

Purpose This User Manual is intended for use with a Milli-Q® Direct Water Purification System. This User Manual is a guide for use during the installation, normal operation and maintenance of a Milli-Q® Direct Water Purification System. It is highly recommended to completely read this manual and to fully comprehend its contents before attempting installation, normal operation or maintenance of the Water Purification System. If this User Manual is not the correct one for your Water Purification System, then please contact Millipore SAS.

Terminology The term "Milli-Q® Direct Water Purification System" is replaced by the terms "Milli-Q® system" or "System" for the remainder of this User Manual unless otherwise noted.

Document FTPF11486 – V3.0, 03/2013

About Millipore SAS

Internet Site Address The Internet site can be used to find addresses, telephone/fax numbers and other information.

Internet Site Address:

www.millipore.com
www.millipore.com/techservice
www.millipore.com/lab_water

Manufacturing Site Millipore SAS
67120 Molsheim
FRANCE

Legal Information

Notice

The information in this document is subject to change without notice and should not be construed as a commitment by Millipore SAS. Millipore SAS assumes no responsibility for any errors that might appear in this document. This manual is believed to be complete and accurate at the time of publication. In no event shall Millipore SAS be liable for incidental or consequential damages in connection with or arising from the use of this manual.

We manufacture and sell water purification systems designed to produce pure or ultra pure water with specific characteristics ($\mu\text{S/cm}$, T, TOC, CFU/ml, Eu/ml) when it leaves the water purification system provided that it's fed with water quality within specifications, and properly maintained as required by the supplier.

We do not warrant these systems for any specific applications. It is up to the end user to determine if the quality of the water produced by our systems matches his expectations, fits with norms/legal requirements and to bear responsibility resulting from the usage of the water.

Product warranty and limitation of liability

The applicable warranty and limitation of liability for the products listed in this publication may be found at <http://www.millipore.com/ec/cp3/terms> within the "Terms and Conditions of Sale" applicable to your purchase transaction.

Trademarks

BioPak, EDS-Pak, Q-Pak, Milli-Q, Millipak, Millipak Express, Progard and Q-POD and are registered trademarks of Merck KGaA, Darmstadt, Germany.

Millipore, the "M" logo, RoClean and VOC-Pak are trademarks of Merck KGaA, Darmstadt, Germany.

All other trademarks are trademarks of their respective owners.

Copyright

© 2013 Millipore SAS, Guyancourt, France. All rights reserved. Printed in France.

This book or parts thereof may not be reproduced in any form without the written permission of the publishers.

The photos are non-contractual.

Recycling



Directive 2002/96 EC: For European users only

The symbol "crossed bin" on a product or its packaging indicates that the product should not be treated like household waste when discarded. Instead the product should be disposed of at a location that handles discarded electric or electronic equipment.

Proper disposal of equipment containing electric or electronic components will help to reduce pollution effects to the environment or to human health. Proper recycling of these products helps in environmental preservation and helps to protect natural resources. For more information about recycling of products containing electric or electronic components, please contact your local recycling representative or organization.

Safety Information

Statement

Your Milli-Q® Direct System should be installed and operated according to the instructions in this manual.

In particular, the hydraulic and electrical specifications should be followed and met. It is important to use this equipment as specified in this manual; using this equipment in a different manner may impair the safety precautions of the Milli-Q® Direct System.

Symbols



This ATTENTION symbol is used to refer to instructions in this manual that need to be done carefully.



These symbols are used to indicate that proper safety equipment has to be used.



Protective glasses and gloves must be worn.



This UV RADIATION sticker is used to refer to a position on the water system Cabinet or inside of it where exposure to UV light is possible.



This DANGER sticker is used to refer to a position on the water system Cabinet or inside of it that could be hazardous.



This ELECTRICAL GROUND sticker is used to refer to a position on the water system Cabinet or inside where an electrical ground connection is made.



This ELECTRICAL DANGER sticker is used to refer to a position on the water system Cabinet or inside where an electrical danger could exist.



IMPORTANT!

Your water system should be installed and operated in a clean and dry area. Please refer to the environment requirements page at the end of this manual.

Your water system is not designed for domestic use.

Table of Contents

Product Information	9
Overview	9
Cabinet	10
Reservoir	15
Consumables	16
Specifications and requirements	17
Installation	20
Overview	20
Alarms generated during installation.....	21
Assembling the POD Unit.....	23
Tubing, cables and power cord.....	24
Installing the Progard® Cartridge.....	27
Installing the Q-Pak® Pack	29
Rinsing the RO Cartridges.....	31
Rinsing the Q-Pak® Pack.....	33
Installing a POD Pak.....	35
Registering UV Lamp timers.....	37
Registering PERFORM RO CL2 CLEANING message timer	39
Registering EXAMINE INLET STRAINER message timer	41
Calibrating the Flow rate	43
Performing a TOC Curve Check.....	45
Software	46
Overview	46
Software Map	47
Standby Mode.....	48
General information.....	48
Description of Standby Menu	49
Manager Menu	52
Description	52
Ready Mode.....	57
General information.....	57
Description of Ready Menu	58
Using the Milli-Q® Direct System	63
Overview	63
Dispensing water.....	64
Viewing water quality	66
Viewing Operation.....	67
Viewing Consumable Status.....	69
Calling Millipore SAS.....	70
Viewing Information	71

Introduction	
Maintenance	73
Overview	73
Maintenance Schedule	74
Replacing the Progard® Cartridge and Vent Filter.....	75
Replacing the Q-Pak® Pack.....	78
Replacing a POD Pak	82
TOC Curve Check.....	85
Sanitising the RO Cartridge(s).....	87
Cleaning the RO Cartridge(s).....	90
Cleaning the Inlet Strainer.....	92
Calibrating the Flow rate.....	95
Alarms	96
Overview.....	96
Alarm Information.....	97
Summary of Alarm messages	101
Alerts	104
Overview.....	104
Alert information.....	105
Summary of Alert messages.....	110
Ordering Information	114
Consumables, Accessories and Systems.....	114

Product Information

Overview

Purpose

This chapter contains topics related to the System.
Some of the more important topics in this chapter are:

- installation requirements,
 - consumable information, and
 - dimensions of various components of the System.
-

Contents

This chapter contains the following topics:

Topic	See Page
Cabinet	10
Reservoir	15
Consumables	16
Specifications and requirements	17

Cabinet

Overview



Item	Description/Name
A	Point Of Delivery (POD)
B	POD Pak
C	Connections for tubings, power cord, level sensor and other cables
D	Q-Pak [®] Pack location
E	Sanitisation Port
F	Main Display
G	Progard [®] Cartridge location


Continued on next page

Cabinet, Continued

Main Display function

The Main Display is used to navigate the System software.

Details of the Main Display



Item	Description
1	LEDs
2	Main LCD
3	Main Keypad



Right

The use of the Right Keypad button is shown below. It is used to move to the next screen.

In this example, the system is changed from STANDBY Mode to READY Mode.

Diagram 1	Action	Diagram 2
<p>Diagram 1 shows the 'STANDBY' screen with the date and time '08 juil. 2009 12:16'. The menu options are 'Menu →' and 'Ready →'. The 'Ready →' option is highlighted in grey and circled in red.</p>	<p>Press .</p>	<p>Diagram 2 shows the 'READY' screen with the date and time '02 juil. 2009 11:45'. The menu options are 'Menu →', 'Tank: Standby →', '0 % Volume →', 'Perm C: 6.0 µS/cm TC', 'MQ Res: - - - MΩcm TC', and 'TOC: - - - ppb'.</p>

Continued on next page

Cabinet, Continued



Left

The use of the Left Keypad button is shown below. It is used to move to the former screen.

Diagram 1	Action	Diagram 2
	Press	



Up

The use of the Up Keypad button is shown below. It is used to scroll up in a menu.

Diagram 1	Action	Diagram 2
	Press	



Down

The use of the Down Keypad button is shown below. It is used to scroll down in a menu.

Diagram 1	Action	Diagram 2
	Press	

Continued on next page

Cabinet, Continued



Validate

The use of the Validate Keypad button is shown below. It is used to confirm a parameter modification.

Diagram 1	Action	Diagram 2
	<p>Press </p>	

READY Mode – water quality values

The READY Mode screen display is explained below.

Diagram	Explanation
	<p>In this example,</p> <ul style="list-style-type: none"> the water filling the tank has a permeate conductivity of 6 $\mu\text{S}/\text{cm}$. the water dispensed from the POD Unit has: <ul style="list-style-type: none"> a resistivity of 18.2 $\text{M}\Omega\cdot\text{cm}$, is temperature compensated (TC) at 25°C, and the TOC value is 4 ppb.
	<p>In this example, there are no Milli-Q® water quality measurements to display. The water quality is only displayed when it is actually measured during water delivery or recirculation.</p>

LEDs

The LEDs are described below.

Item	Description
Green LED	System is operating within specifications.
Yellow LED	An Alert is present.
Red LED	An Alarm is present.


NOTE:

If an Alarm and an Alert are present at the same time, then only the red LED is lit. The red and yellow LEDs are never lit at the same time.

Continued on next page

Cabinet, Continued

Port and cables The port and cable connections are explained below.



Item	Description	Item	Description
1	RO Reject Port	A	Power Entry connection (100 – 240 V)
2	Feed water Port	B	Accessories connection (maximum 24 VDC)
3	From Reservoir Port	C	Termination Plug Connection
4	To Reservoir Port	D	Level Sensor Connection (maximum 5 VDC)
		E	Ethernet connection (maximum 5 VDC)

Reservoir

Information


Millipore SAS recommends using a Reservoir having the following catalogue number:

Size	Catalogue Numbers
30 Litre	TANKPE030
60 Litre	TANKPE060
100 Litre	TANKPE100

Consumables

Flow diagram

The water flow through a System is shown here in a flow diagram. The various consumables are described below.



Item	Description
A	Progard® Cartridge
B	UV 254 nm Lamp (Optional)
C	UV 185 nm Lamp
D	Q-Pak® Pack
E	POD Pak

Progard® Cartridge

The Progard® Cartridge protects the RO Cartridge in order to increase its lifetime. It prevents mineral scaling, damage due to particulate and chlorine oxidation of the RO Cartridge(s).

UV 254 nm Lamp

The UV 254 nm Lamp is optional. It emits light at 254 nm. It is used to kill bacteria.

UV 185 nm Lamp

The dual wavelength UV 185 nm Lamp emits light at 185 nm and at 254 nm. It kills bacteria and reduces the level of organic molecules in the water.

Q-Pak® Pack

The Q-Pak® Pack removes trace levels of ions and organic molecules.

POD Pak

The POD Pak is the final water purification device. It is attached to the Point of Delivery outlet. The POD Pak provides additional quality and insurance that trace contaminants related to specific applications are removed just before ultrapure water is delivered.

Specifications and requirements

Milli-Q® Water quality

The water delivered from a POD Unit has the following characteristics.

Parameter	Specification	Units
Resistivity	18.2	MΩ.cm @25°C
TOC	≤ 5	ppb
Particulates > 0.22 µm**	< 1	Particulates/mL
Bacteria**	< 0.1	cfu/mL
Pyrogens*	< 0.001	Eu/mL
RNases*	< 0.01	ng/mL
DNases*	< 4	pg/µL
Flow Rate**	0.05 – 1.5	L/min

(*) With BioPak® Final Filter

(**) With Millipak® or BioPak® Final Filter

NOTE:

These specifications are valid if feed water within specification and if correct maintenance is performed on the system. Some specifications may not be achieved at start-up.

Weight

The various weights are found in the table below.

System	Operating Weight (kg)	Dry Weight (kg)	Shipping Weight (kg)
Milli-Q® Direct 8	27	20	24
Milli-Q® Direct 16	28	21	25

Electrical


The electrical specifications and data are found in the table below.

Parameter	Value
Voltage	100-230 VAC ±10%
Frequency	50-60 Hz ±10%
Main Fuse	<ul style="list-style-type: none"> • 3.15 Amp Fast Acting; 5 mm x 20 mm; 250 V safety voltage. • The fuse should be serviced by a qualified Millipore SAS Service Representative.
Power Used	145 VA
Power Cord Length	2.5 metres
Electrical Ground	Earth Grounded
Power Cord use	<ul style="list-style-type: none"> • The System is powered on and off by removing the power cord from the wall outlet. • The power cord should be plugged into a wall outlet that is accessible.

Continued on next page

Specifications and requirements, Continued

Dimensions



Materials of construction

Please contact Millipore SAS for a list of the Materials of Construction.

Continued on next page

Specifications and requirements, Continued

Feed water The Feed water requirements are listed here.

Parameter	Value
Type of Feed water	Potable tap water
Conductivity	< 2000 $\mu\text{S}/\text{cm}$
Pressure	1 bar < P < 6 bar
Temperature	5°C < T < 35°C
Dissolved CO ₂	< 30 ppm
Free Chlorine	< 3 ppm
Fouling Index	< 12
pH	4 < pH < 10

Environmental The Environmental requirements are listed here.

Parameter	Value
Altitude	< 3000 metres
Ambient operating temperature	4 – 40°C
Ambient storage temperature	4 – 40°C
Installation Category	II
Location	The System is intended for indoor use only.
Pollution Degree	2
Relative humidity during storage and operation	Maximum relative humidity 80% for temperatures up to 31°C decreasing linearly to 50% relative humidity at 40°C.

Noise Level The noise level is < 50 dB at a distance of 1 metre.

Consumables The minimum consumables required for installation are listed here. Note that these items are not shipped with the System and must be ordered separately:

- Progard® Cartridge,
- Q-Pak® Pack, and
- POD Pak.

Reservoir location The Reservoir must be located relative to the Water System:

- $0 \leq y \leq 2$ metres, where y = vertical distance, and
- $0 \leq x \leq 3$ metres, where x = horizontal distance.

Installation

Overview

Purpose This chapter explains how to install the System.

Contents This chapter contains the following topics:

Topic	See Page
Alarms generated during installation	21
Assembling the POD Unit	23
Tubing, cables and power cord	24
Installing the Progard® Cartridge	27
Installing the Q-Pak® Pack	29
Rinsing the RO Cartridges	31
Rinsing the Q-Pak® Pack	33
Installing a POD Pak	35
Registering UV Lamp timers	37
Registering PERFORM RO CL2 CLEANING message timer	39
Registering EXAMINE INLET STRAINER message timer	41
Calibrating the Flow rate	43
Performing a TOC Curve Check	45

Summary list The steps shown below outline the sequence and major actions of a System installation. Please refer to this list throughout the installation.

Step	Action
1	Put POD Arm onto POD Mast
2	Put Point Of Delivery onto POD Arm
3	Install tubing, termination plug and power cord
4	Power on the System, check date and time
5	Install the Q-Pak® Pack
6	Install, and flush the Progard® Cartridge
7	Flush and rinse the RO Cartridge(s)
8	Fill the Reservoir
9	Flush and rinse the Q-Pak® Pack
10	Install and Register the POD Pak
11	Register the UV Lamp timers
12	Register the PERFORM RO CL2 CLEANING message timer
13	Register the EXAMINE INLET STRAINER message timer
14	Calibrate the Product Water flow rate
15	Perform a TOC Curve Check

Alarms generated during installation

Overview

- During the installation of a Milli-Q® System, certain Alarm messages are generated.
 - This occurs because:
 - the Reservoir is empty,
 - there is air in the tubings and in the Progard® Cartridge,
 - the Progard® Cartridge is not installed, and
 - the Q-Pak® Pack is not installed.
 - These alarms are explained here. For more information about Alarm messages, see the chapter titled 'Alarms'.
-

TANK EMPTY message

- This alarm occurs because the Reservoir is empty during most of the installation.
 - This alarm goes away when the Reservoir is partially full.
 - To cancel the text display of this alarm message, follow the instructions on the LCD.
-

PROGARD CARTRIDGE OUT message

- This alarm occurs because the Progard® Cartridge is not installed.
 - This alarm goes away when the Progard® Cartridge is detected by the Milli-Q® System.
 - To cancel the text display of this alarm message, follow the instructions on the LCD.
-

Q-PAK PACK OUT message

- This alarm occurs because the Q-Pak® Pack is not installed.
 - This alarm goes away when the Q-Pak® Pack is detected by the System.
 - To cancel the text display of this alarm message, follow the instructions on the LCD.
-

MILLI-Q RES < SP, REPLACE Q-PAK message

- This alarm occurs because the Q-Pak® Pack is not fully rinsed out or there is air in the tubing near the resistivity sensor.
 - This alarm goes away when a few litres of water are dispensed from the POD Unit.
 - To cancel the text display of this alarm message, follow the instructions on the LCD.
-

LOW FEED WATER PRESSURE message

- This alarm occurs because there is air in the tubings and in the new Progard® cartridge.
 - When the air is gone and replaced with water, this alarm does not occur anymore during installation.
 - To cancel the text display of this alarm message, follow the instructions on the LCD.
-

Continued on next page

Alarms generated during installation, Continued


MILLI-Q TOC >
SP

- This alarm occurs because the TOC indicator algorithm needs data after the Q-Pak® pack is fully rinsed out.
 - To cancel the text display of this alarm message, follow the instructions on the LCD.
-

Assembling the POD Unit

Separating POD Arm and Point Of Delivery

Separate the POD Arm and the Point Of Delivery by cutting and removing the tape that holds them together.



Placing the POD Arm


Place the POD and POD Arm onto the POD Mast as shown below.



Tubing, cables and power cord

Summary

Item	Description
1	RO Reject Water tubing. Goes to a drain.
2	Feed water supply tubing to Milli-Q® Direct system.
3	Tubing connected here comes from the bottom of the Reservoir. See the next section.
4	Tubing connected here goes to the bottom of the Reservoir. See the next section.
A	Power cord
B	Accessories cable
C	Termination Plug
D	Level Sensor from Reservoir
E	Ethernet cable




Continued on next page

Tubing, cables and power cord, Continued


Feed water tubing to pipe

- Install the Inlet Strainer as shown here.
- Connect one end of the feed water tubing to the Inlet Strainer.



Reservoir connections

The tubings from the Water System Ports 3 and 4 are connected to the Reservoir as shown here.



NOTE:

The valve where the tubing from Port 3 is connected must be opened.

Powering the system

- Open the feed water source.
 - Plug the power cord into the Water System.
 - Plug the power cord into a source of electrical power.
 - The Main LCD shows a series of start-up screens.
-

Continued on next page

Tubing, cables and power cord, Continued

Alarm messages Because the System is starting with an empty tank, without a Progard® Cartridge or a Q-Pak® Pack installed, there are alarm messages displayed.

These alarms are:

- TANK EMPTY,
 - Q-PAK® PACK OUT, and
 - PROGARD CARTRIDGE OUT.
-

Cancel Alarms When an Alarm message is displayed, follow the instructions on the screen to cancel the text display of the Alarm.





Check the date

- When the Alarm messages are cancelled, check that the displayed date is correct.
 - If necessary, go to the Manager Menu Software and correct the date and time. See the Software Map in the beginning of the Software Chapter for more information.
 - Do not install a Progard® Cartridge or a Q-Pak® Pack until the displayed date is correct.
-

Installing the Progard® Cartridge

Procedure

Follow the steps below to install a new Progard® Cartridge.

Step	Action	Diagram
1	<p>Start in STANDBY Mode.</p> <p>NOTE: The PROGARD CARTRIDGE OUT Alarm message is not shown at this time. By following the instructions earlier in this manual, the alarm was cancelled.</p>	
2	<ul style="list-style-type: none"> • Open the right door of the System Cabinet. • Remove the 2 protective caps located on the ports inside. 	
3	<ul style="list-style-type: none"> • Remove the covers on the 2 ports of the Progard® Cartridge. • Wet the O-rings with water. 	
4	<ul style="list-style-type: none"> • Install the Progard® Cartridge until it is fully seated. • Close the right door. 	
5	<p>One minute later, the Main LCD shows that a new Progard® Cartridge is installed.</p>	

Continued on next page

Installing the Progard® Cartridge, Continued


Procedure
(continued)

Step	Action	Diagram
6	Press  .	
7	When the Progard® Cartridge flush has finished, the Water System goes to READY Mode.	

Installing the Q-Pak® Pack

Procedure

Follow the steps below to install a new Q-Pak® Pack.

Step	Action	Diagram
1	Start in STANDBY Mode.	
2	<ul style="list-style-type: none"> • Open the left door of the System. • Remove the 2 protective caps located on the ports inside. 	
3	<ul style="list-style-type: none"> • Remove the covers on the 2 ports of the Q-Pak® Pack. • Make sure the rubber O-rings are firmly in place. • Wet the O-rings with water. 	
4	Push the top of the Q-Pak® Pack into the ports on the System.	

Continued on next page

Installing the Q-Pak® Pack, Continued

Procedure
(continued)

Step	Action	Diagram
5	Push the bottom of the Q-Pak® Pack inwards.	
6	Push the pack locking handle down. Close the left door.	
7	One minute later, the Main LCD shows that a new Q-Pak® Pack is installed.	
8	Press  .	

Rinsing the RO Cartridges



Rinse the RO Cartridges

The RO Cartridges must be flushed and rinsed when the Milli-Q® System is installed. Failure to do this results in poor water quality.

Procedure

Follow the steps below to flush and rinse the RO Cartridge(s).

Step	Action	Diagram
1	Start in STANDBY Mode.	
2	<ul style="list-style-type: none"> Select Menu. Press 	
3	<ul style="list-style-type: none"> Select Maintenance. Press 	
4	<ul style="list-style-type: none"> Select Install new RO. Press 	

Continued on next page

Rinsing the RO Cartridges, Continued

Procedure
(continued)

Step	Action	Diagram
5	Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>INSTALL NEW RO</p> <p>This procedure should be performed by a Millipore trained service engineer. Press → to continue or ← to exit.</p> </div>
6	Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>INSTALL NEW RO</p> <p>The Millipore trained service engineer confirms RO cartridge installation by pressing ✓. A 15 minute RO Flush Followed by a 225 minute RO rinse will start. Press ← to exit.</p> </div>
7	Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>INSTALL NEW RO</p> <p>RO Flush in progress. Remaining Time : 15 min.</p> </div>
8	After 15 minutes, the LCD looks like this.	<div style="border: 1px solid black; padding: 5px;"> <p>INSTALL NEW RO</p> <p>RO Rinse in progress. Remaining Time : 225 min.</p> </div>
9	When the 225 minute RO rinse is finished, the Milli-Q® System returns to READY Mode. The Reservoir is now being filled.	<div style="border: 1px solid black; padding: 5px;"> <p>READY</p> <p>02 jul. 2009 11:48</p> <p style="text-align: right;">Menu →</p> <p>Tank: Standby →</p> <p>0 % Volume →</p> <p>Perm C: 6.0 µS/cm TC</p> <p>MQ Res: - - - MΩ·cm TC</p> <p>TOC: - - - ppb</p> </div>

Rinsing the Q-Pak® Pack



Have enough water!

There has to be enough water in the Reservoir in order to rinse the Q-Pak® Pack.

Reservoir	Minimum
30 Litre	100% Tank Level
60 Litre	> 40% Tank Level
100 Litre	> 30% Tank Level

If there is not enough water in the Reservoir, the TANK EMPTY Alarm is shown. Additionally, air can enter the tubings and can temporarily affect other sensors.

Procedure

Follow the steps below to rinse the Q-Pak® Pack.

Step	Action	Diagram
1	<ul style="list-style-type: none"> • Locate the clear tubing and the barbed fitting from the System Accessories Bag. • Screw the barbed fitting onto the POD Unit. • Push one end of the clear tubing onto the end of the barbed fitting. • Place the other end of the clear tubing into a sink. <p>NOTE: Do not use any white tape on the threads of the barbed fitting. An O-ring located inside the POD Dispenser ensures water tightness.</p>	
2	Place the System into READY Mode.	
3	Push the POD Plunger all the way down and then release it. In a few minutes, water should come out of the POD Unit.	

Continued on next page

Rinsing the Q-Pak® Pack, Continued

Procedure (continued)

Step	Action	Diagram
4	Dispense water for at least 10 minutes.	 <p>READY 02 juil. 2009 13:55 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 MΩcm TC TOC: 4 ppb</p>
5	Push the POD Plunger all the way down and then release it to stop dispensing water. Leave the System in READY Mode.	 <p>READY 02 juil. 2009 13:55 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 MΩcm TC TOC: 4 ppb</p>

Installing a POD Pak

Overview

The installation of a POD Pak involves 2 steps. These are:


- placing and flushing the POD Pak onto the POD Unit, and
- registering the installation of a specific POD Pak.

Placing and flushing

Follow the instructions delivered with the POD Pak.

Registering

Follow the steps below to register the installation of the POD Pak.

Step	Action	Diagram
1	Start in STANDBY Mode.	
2	Select Menu. Press  .	
3	Select Maintenance. Press  .	
4	Scroll down to Install POD Pak. Select it.	

Continued on next page

Installing a POD Pak, Continued

Registering (continued)

Step	Action	Diagram
5	Press  .	
6	Press  .	
7	In this example, you choose Millipak®. Press  .	
8	Press  .	
9	Press  .	
10	Press 3 times on  .	

Registering UV Lamp timers

Introduction

The timer used for each UV Lamp must be reset when the System is installed. If this is not done, then the message indicating that a Lamp replacement is needed is shown too early.

The UV Lamp timers need to be reset for:


- the UV 185 nm Lamp, and
- the UV 254 nm Lamp.

NOTE:

Before doing this, make sure that the date and time have been checked for accuracy.

Procedure

This procedure shows how to reset the timer used for the UV 185 nm Lamp.

Step	Action	Diagram
1	Place the System in STANDBY Mode.	
2	Select Menu. Press  .	
3	Select Maintenance. Press  .	
4	Select Install UV 185 nm Lamp. Press  .	

Continued on next page

Registering UV Lamp timers, Continued

Procedure
(continued)

Step	Action	Diagram
5	Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>INSTALL UV 185 LAMP</p> <p>This procedure should be performed by a Millipore trained service engineer. Press → to continue or ← to exit.</p> </div>
6	Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>INSTALL UV 185 LAMP</p> <p>The Millipore trained service engineer confirms the UV 185 nm Lamp installation by pressing ✓. Press ← to exit.</p> </div>
7	Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>INSTALL UV 185 LAMP</p> <p>UV 185 nm Lamp installation is registered. Next maintenance in 730 days. Press ← to exit.</p> </div>
8	Press 3 times on  .	<div style="border: 1px solid black; padding: 5px;"> <p>STANDBY</p> <p>02 juil. 2009 13:57</p> <p style="text-align: right;">Menu →</p> <p style="text-align: right;">Ready →</p> </div>

**Reset UV
254nm Lamp
timer**

After resetting the timer for the UV 185 nm Lamp timer, reset the UV Lamp timer for the UV 254 nm Lamp.

Registering PERFORM RO CL2 CLEANING message timer

Introduction



- The timer used to perform RO CL2 cleaning must be reset when the Water System is installed.
- If this is not done, then the message indicating that the message PERFORM RO CL2 CLEANING is shown too early.

Note

This is only done once, at installation. In the future, this timer is reset automatically after an RO CL2 cleaning is performed.

Procedure

This procedure shows how to reset the timer used for the message PERFORM RO CL2 CLEANING.

Step	Action	Diagram
1	Start in STANDBY Mode.	
2	<ul style="list-style-type: none"> • Select Menu. • Press . 	
3	<ul style="list-style-type: none"> • Select Maintenance. • Press . 	
4	<ul style="list-style-type: none"> • Select Reset RO CL2 Cleaning. • Press . 	

Continued on next page

Registering PERFORM RO CL2 CLEANING message timer, Continued

Procedure
(continued)

Step	Action	Diagram
5	Press  .	
6	Press 3 times on  .	


Registering EXAMINE INLET STRAINER message timer

Introduction

- The timer used for cleaning the Inlet Strainer must be reset when the Water System is installed.
- If this is not done, then the message indicating that the message EXAMINE INLET STRAINER is shown too early.

Procedure

This procedure shows how to reset the timer used for the message EXAMINE INLET STRAINER.

Step	Action	Diagram
1	Start in STANDBY Mode.	
2	<ul style="list-style-type: none"> • Select Menu. • Press . 	
3	<ul style="list-style-type: none"> • Select Maintenance. • Press . 	
4	<ul style="list-style-type: none"> • Select Clean Strainer. • Press . 	

Continued on next page

Registering EXAMINE INLET STRAINER message timer, Continued




Procedure
(continued)

Step	Action	Diagram
5	Press  .	<div data-bbox="970 342 1230 607" style="border: 1px solid black; padding: 5px;"> <p>CLEAN STRAINER</p> <p>See Maintenance Chapter in the User Manual For more information.</p> <p>Press ✓ after cleaning or ← to exit.</p> </div>
6	Press  .	<div data-bbox="970 622 1230 887" style="border: 1px solid black; padding: 5px;"> <p>CLEAN STRAINER</p> <p>The strainer cleaning date is registered. Next maintenance in 365 days.</p> <p>Press ← to exit.</p> </div>
7	Press 3 times on  .	<div data-bbox="970 902 1230 1167" style="border: 1px solid black; padding: 5px;"> <p>STANDBY</p> <p>02 juil. 2009 13:57</p> <p style="text-align: right;">Menu →</p> <p style="text-align: right;">Ready →</p> </div>

Calibrating the Flow rate

Introduction The Milli-Q® Water flow rate should be calibrated when the System is installed. A 1 Litre graduated cylinder is needed.

Procedure Follow the steps below to perform a Flow Calibration.

Step	Action	Diagram
1	Go to STANDBY Mode.	
2	Select Menu. Press  .	
3	Enter the Manager Menu. See the Software Chapter to learn how to enter the Manager Menu.	
4	Select Setup. Press  .	
5	Select Flow Calibration. Press  .	

Continued on next page

Calibrating the Flow rate, Continued

Procedure
(continued)

Step	Action	Diagram
6	Place a 1 L Graduated Cylinder under the POD Unit. Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>FLOW CALIBRATION</p> <p>Press ✓ or press  on the Q-POD keypad if you have installed one to start water delivery.</p> <p>After the water dispensing is complete, measure the collected volume.</p> </div>
7	Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>FLOW CALIBRATION</p> <p>The system is now delivering water.</p> <p>Task Completion: XX %</p> </div>
8	Water dispenses automatically from the POD Unit. Wait until it stops dispensing water.	<div style="border: 1px solid black; padding: 5px;"> <p>FLOW CALIBRATION</p> <p>Volume : 900 mL</p> <p>Use + and - keys to register the value of the collected volume. Press ✓ to confirm and exit.</p> </div>
9	Measure the amount of water (in ml) that was dispensed. Suppose 870 ml was collected. Input this using the Keypad.	<div style="border: 1px solid black; padding: 5px;"> <p>FLOW CALIBRATION</p> <p>Volume : 870 mL</p> <p>Use + and - keys to register the value of the collected volume. Press ✓ to confirm and exit.</p> </div>
10	Perform the flow calibration again to improve accuracy. Press  .	<div style="border: 1px solid black; padding: 5px;"> <p>SETUP</p> <p>Install Date →</p> <p>Buzzer →</p> <p>MQ Recirc Mode →</p> <p>POD Flow Stop →</p> <p>Temp Comp Mode →</p> <p>Flow Calibration →</p> <p>UV 254 nm Activation →</p> </div>
11	Press 3 times on  .	<div style="border: 1px solid black; padding: 5px;"> <p>STANDBY</p> <p>02 juil. 2009 14:27</p> <p>Menu →</p> <p>Ready →</p> </div>

Performing a TOC Curve Check

Introduction

The indication of TOC values is performed with information supplied with a TOC Curve Check. In order to update this information, perform a TOC Curve Check by following the steps below.

Procedure

Follow the steps below to perform a TOC Curve Check.

Step	Action	Diagram
1	Go to READY Mode.	
2	Select Menu. Press	
3	Select TOC Curve Check. Press	
4	Press	
5	After approximately 10 minutes, the System returns to READY Mode.	

Software

Overview

Introduction


The purpose of this chapter is to explain the various software used in the System.

Contents

This chapter contains the following topics:

Topic	See Page
Software Map	47
Standby Mode	48
Manager Menu	52
Ready Mode	57

Software Map



Standby Mode

General information

Purpose

STANDBY mode is used primarily for:

- maintenance actions, and
 - going to the Manager Menu.
-

Display




READY Mode from STANDBY Mode

Diagram 1	Action	Diagram 2
	<p>Press .</p>	

Description of Standby Menu

Maintenance The Maintenance Menu is described below.



Item	Description
Install Pretreatment	Used to reset Alert message REPLACE EXTERNAL PRE-TREATMENT.
Clean Strainer	Used to reset Alert message EXAMINE INLET STRAINER.
Install Progard®	Used to see general information about the Progard® Cartridge exchange.
Install new RO	Used to start a flush and rinse of a new RO Cartridge.
Install UV 254 Lamp	Used to reset Alert message REPLACE 254 NM LAMP.
Install UV 185 Lamp	Used to reset Alert message REPLACE 185 NM LAMP.
Install Q-Pak®	Used to see general information about the Q-Pak® pack exchange.
Install POD Pak	Used to reset Alert message REPLACE POD PAK.
Install ASM UV	Used to reset Alert message REPLACE ASM UV LAMP
Reset RO CL2 Cleaning	Used to reset Alert message PERFORM RO CL2 CLEANING at installation.

Continued on next page

Description of Standby Menu, Continued

Sanitise/clean

Diagram 1	Diagram 2															
<table border="1"> <tr><td>STANDBY MENU</td></tr> <tr><td>Maintenance →</td></tr> <tr><td>Sanitise/Clean →</td></tr> <tr><td>Suitability Tests →</td></tr> <tr><td>Language →</td></tr> <tr><td>Manager Menu →</td></tr> <tr><td> </td></tr> <tr><td> </td></tr> </table>	STANDBY MENU	Maintenance →	Sanitise/Clean →	Suitability Tests →	Language →	Manager Menu →			<table border="1"> <tr><td>SANITISE / CLEAN</td></tr> <tr><td>RO CL2 Cleaning →</td></tr> <tr><td>RO pH Cleaning →</td></tr> <tr><td>RO Cleaning →</td></tr> <tr><td>System Cleaning →</td></tr> <tr><td> </td></tr> <tr><td> </td></tr> </table>	SANITISE / CLEAN	RO CL2 Cleaning →	RO pH Cleaning →	RO Cleaning →	System Cleaning →		
STANDBY MENU																
Maintenance →																
Sanitise/Clean →																
Suitability Tests →																
Language →																
Manager Menu →																
SANITISE / CLEAN																
RO CL2 Cleaning →																
RO pH Cleaning →																
RO Cleaning →																
System Cleaning →																

Item	Description
RO CL2 Cleaning	Used to sanitise the RO Cartridge(s).
RO pH Cleaning	Used to clean the RO Cartridge(s).
System Cleaning	Contact Millipore SAS for more information.

Suitability Tests

Diagram 1	Diagram 2													
<table border="1"> <tr><td>STANDBY MENU</td></tr> <tr><td>Maintenance →</td></tr> <tr><td>Sanitise/Clean →</td></tr> <tr><td>Suitability Tests →</td></tr> <tr><td>Language →</td></tr> <tr><td>Manager Menu →</td></tr> <tr><td> </td></tr> <tr><td> </td></tr> </table>	STANDBY MENU	Maintenance →	Sanitise/Clean →	Suitability Tests →	Language →	Manager Menu →			<table border="1"> <tr><td>SUITABILITY TESTS</td></tr> <tr><td>Res Suitability Test →</td></tr> <tr><td>Temp Suitability Test →</td></tr> <tr><td> </td></tr> <tr><td> </td></tr> </table>	SUITABILITY TESTS	Res Suitability Test →	Temp Suitability Test →		
STANDBY MENU														
Maintenance →														
Sanitise/Clean →														
Suitability Tests →														
Language →														
Manager Menu →														
SUITABILITY TESTS														
Res Suitability Test →														
Temp Suitability Test →														

Item	Description
Res Suitability Test	Contact Millipore SAS for more information.
Temp Suitability Test	

Continued on next page

Description of Standby Menu, Continued

Language

Diagram 1	Diagram 2																	
<table border="1"> <tr><td>STANDBY MENU</td></tr> <tr><td>Maintenance →</td></tr> <tr><td>Sanitise/Clean →</td></tr> <tr><td>Suitability Tests →</td></tr> <tr style="background-color: #cccccc;"><td>Language →</td></tr> <tr><td>Manager Menu →</td></tr> <tr><td> </td></tr> <tr><td> </td></tr> </table>	STANDBY MENU	Maintenance →	Sanitise/Clean →	Suitability Tests →	Language →	Manager Menu →			<table border="1"> <tr><td>LANGUAGE</td></tr> <tr><td>Chinese</td></tr> <tr style="background-color: #cccccc;"><td>English ✓</td></tr> <tr><td>French</td></tr> <tr><td>German</td></tr> <tr><td>Italian</td></tr> <tr><td>Japanese</td></tr> <tr><td>Portuguese</td></tr> <tr><td> </td></tr> </table>	LANGUAGE	Chinese	English ✓	French	German	Italian	Japanese	Portuguese	
STANDBY MENU																		
Maintenance →																		
Sanitise/Clean →																		
Suitability Tests →																		
Language →																		
Manager Menu →																		
LANGUAGE																		
Chinese																		
English ✓																		
French																		
German																		
Italian																		
Japanese																		
Portuguese																		

Item	Description
Language	Change the displayed language.

Manager Menu See the next section for information about the Manager Menu.


Manager Menu

Description

How to enter


- See the Software Map at the beginning of this chapter. The map shows how to enter the Manager Menu.
- To enter the Manager Menu, it is necessary to input a Login and a Password.
- The Software Map indicates how to input a Login and a Password.

Change ID and Password



Item	Description
Change ID & Password	<ul style="list-style-type: none"> • Change the Login and Password used to enter the Manager Menu. • Use 4 characters for the Login and the Password.

Date and Time



Item	Description
Date and Time	Adjust your local date and time.

Continued on next page

Description, Continued

Set Points

Diagram 1	Diagram 2																									
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Change ID and Password →</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> </table>	MANAGER MENU	Change ID and Password →	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	<table border="1"> <tr><td>SET POINTS</td></tr> <tr><td>Pretreatment →</td></tr> <tr><td>Strainer Frequency →</td></tr> <tr><td>Tap Feed Cond →</td></tr> <tr><td>RO Rejection →</td></tr> <tr><td>RO CL2 Cleaning →</td></tr> <tr><td>Permeate Cond →</td></tr> <tr><td>Tank Refill →</td></tr> </table>	SET POINTS	Pretreatment →	Strainer Frequency →	Tap Feed Cond →	RO Rejection →	RO CL2 Cleaning →	Permeate Cond →	Tank Refill →	<table border="1"> <tr><td>SET POINTS</td></tr> <tr><td>Milli-Q Product Res →</td></tr> <tr><td>Milli-Q Product TOC →</td></tr> <tr><td>Millipak →</td></tr> <tr><td>BioPak →</td></tr> <tr><td>EDS-Pak →</td></tr> <tr><td>Pod Pak A →</td></tr> <tr><td>Pod Pak B →</td></tr> </table>	SET POINTS	Milli-Q Product Res →	Milli-Q Product TOC →	Millipak →	BioPak →	EDS-Pak →	Pod Pak A →	Pod Pak B →
MANAGER MENU																										
Change ID and Password →																										
Date and Time →																										
Set Points →																										
Units →																										
Setup →																										
User Parameters →																										
History →																										
SET POINTS																										
Pretreatment →																										
Strainer Frequency →																										
Tap Feed Cond →																										
RO Rejection →																										
RO CL2 Cleaning →																										
Permeate Cond →																										
Tank Refill →																										
SET POINTS																										
Milli-Q Product Res →																										
Milli-Q Product TOC →																										
Millipak →																										
BioPak →																										
EDS-Pak →																										
Pod Pak A →																										
Pod Pak B →																										

Item	Description
Pretreatment	Change set point for controlling the frequency of the message REPLACE EXTERNAL PRE-TREATMENT.
Strainer Frequency	Change set points for controlling the frequency of the message EXAMINE INLET STRAINER.
Tap Feed Cond	Change set point controlling the message TAP FEED CONDUCTIVITY > SP.
RO Rejection	Change set point controlling the message RO REJECTION < SP.
RO CL2 Cleaning	Change set point for controlling the frequency of the message PERFORM RO CL2 CLEANING
Permeate Cond	Change set point controlling the message PERMEATE C > SP.
Tank Refill	Change set point controlling the tank level where the Milli-Q® System starts to refill the tank.
Milli-Q® Product Res	Change set point controlling the message MILLI-Q RES < SP, REPLACE Q-PAK.
Milli-Q® Product TOC	Change set point controlling the message MILLI-Q TOC > SP.
Millipak®	Change set point controlling the message REPLACE POD PAK IN XX DAYS (where $1 \leq XX \leq 15$).
BioPak®, EDS-Pak®, POD Pak	See above.

Continued on next page

Description, Continued

Units

Diagram 1	Diagram 2															
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Change ID and Password →</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> <tr><td> </td></tr> </table>	MANAGER MENU	Change ID and Password →	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →		<table border="1"> <tr><td>UNITS</td></tr> <tr><td>Pressure →</td></tr> <tr><td>Milli-Q Product →</td></tr> <tr><td>Tank Volume →</td></tr> <tr><td> </td></tr> <tr><td> </td></tr> </table>	UNITS	Pressure →	Milli-Q Product →	Tank Volume →		
MANAGER MENU																
Change ID and Password →																
Date and Time →																
Set Points →																
Units →																
Setup →																
User Parameters →																
History →																
UNITS																
Pressure →																
Milli-Q Product →																
Tank Volume →																

Item	Description
Pressure	<ul style="list-style-type: none"> • Change the displayed units of pressure. • Choices are bar, psi and KPa.
Milli-Q® Product	<ul style="list-style-type: none"> • Change the displayed units of Milli-Q® Product Water quality. • Choices are MΩ.cm or μS/cm.
Tank Volume	<ul style="list-style-type: none"> • Change the displayed units of Tank Volume. • Choices are % full, Litres or US Gallons.

Continued on next page

Description, Continued

Setup

Diagram 1	Diagram 2																								
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Change ID and Password →</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> </table>	MANAGER MENU	Change ID and Password →	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	<table border="1"> <tr><td>SETUP</td></tr> <tr><td>Install Date →</td></tr> <tr><td>Buzzer →</td></tr> <tr><td>MQ Recirc Mode →</td></tr> <tr><td>POD Flow Stop →</td></tr> <tr><td>Temp Comp Mode →</td></tr> <tr><td>Flow Calibration →</td></tr> <tr><td>UV 254 nm Activation →</td></tr> </table>	SETUP	Install Date →	Buzzer →	MQ Recirc Mode →	POD Flow Stop →	Temp Comp Mode →	Flow Calibration →	UV 254 nm Activation →	<table border="1"> <tr><td>SETUP</td></tr> <tr><td>Flow Calibration →</td></tr> <tr><td>UV 254 nm Activation →</td></tr> <tr><td>UV 185 nm Activation →</td></tr> <tr><td>ASM UV Lamp Schedule →</td></tr> <tr><td>Network Settings →</td></tr> <tr><td>TOC Curve Check Time →</td></tr> </table>	SETUP	Flow Calibration →	UV 254 nm Activation →	UV 185 nm Activation →	ASM UV Lamp Schedule →	Network Settings →	TOC Curve Check Time →
MANAGER MENU																									
Change ID and Password →																									
Date and Time →																									
Set Points →																									
Units →																									
Setup →																									
User Parameters →																									
History →																									
SETUP																									
Install Date →																									
Buzzer →																									
MQ Recirc Mode →																									
POD Flow Stop →																									
Temp Comp Mode →																									
Flow Calibration →																									
UV 254 nm Activation →																									
SETUP																									
Flow Calibration →																									
UV 254 nm Activation →																									
UV 185 nm Activation →																									
ASM UV Lamp Schedule →																									
Network Settings →																									
TOC Curve Check Time →																									

Item	Description
Install Date	Change the installation date.
Buzzer	Change the setting for the Buzzer.
MQ Recirc Mode	Change the amount of time that the System automatically recirculates every hour in READY Mode. NOTE If set to 60 minutes, the daily TOC Curve Check will not be performed.
POD Flow Stop	Change the amount of time that the POD Unit dispenses continuously before it automatically stops.
Temp Comp Mode	Change the Temperature Compensation Mode.
Flow Calibration	Used for performing a flow calibration.
UV 254 nm Activation	Used to activate or deactivate the UV 254 nm Lamp.
UV 185 nm Activation	Used to activate or deactivate the UV 185 nm Lamp.
ASM UV Lamp Schedule	<ul style="list-style-type: none"> Used to change the times when the ASM (Automatic Sanitisation Module) turns on. See the ASM User Manual for more information.
Network Settings	<ul style="list-style-type: none"> Change Network settings. Contact Millipore SAS for more information.
TOC Curve Check Time	Change the time when a TOC Curve Check is automatically performed each day.

Continued on next page

Description, Continued

User Parameters The User Parameters are seen when a History Report is printed out.

Diagram 1	Diagram 2																
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Change ID and Password →</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> </table>	MANAGER MENU	Change ID and Password →	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	<table border="1"> <tr><td>USER PARAMETERS</td></tr> <tr><td>Company Name →</td></tr> <tr><td>Department Name →</td></tr> <tr><td>Address →</td></tr> <tr><td>Postal Code →</td></tr> <tr><td>City →</td></tr> <tr><td>Country →</td></tr> <tr><td>Email →</td></tr> </table>	USER PARAMETERS	Company Name →	Department Name →	Address →	Postal Code →	City →	Country →	Email →
MANAGER MENU																	
Change ID and Password →																	
Date and Time →																	
Set Points →																	
Units →																	
Setup →																	
User Parameters →																	
History →																	
USER PARAMETERS																	
Company Name →																	
Department Name →																	
Address →																	
Postal Code →																	
City →																	
Country →																	
Email →																	

Item	Description
Company Name	Change the item
Department Name	
Address	
Postal Code	
City	
Country	
Email	

History Summary

Diagram 1	Diagram 2														
<table border="1"> <tr><td>MANAGER MENU</td></tr> <tr><td>Date and Time →</td></tr> <tr><td>Set Points →</td></tr> <tr><td>Units →</td></tr> <tr><td>Setup →</td></tr> <tr><td>User Parameters →</td></tr> <tr><td>History →</td></tr> <tr><td>Lab closed →</td></tr> </table>	MANAGER MENU	Date and Time →	Set Points →	Units →	Setup →	User Parameters →	History →	Lab closed →	<table border="1"> <tr><td>HISTORY</td></tr> <tr><td>History Summary →</td></tr> <tr><td>Print System History →</td></tr> <tr><td>Print RO History →</td></tr> <tr><td>Print Milli-Q History →</td></tr> <tr><td>Print Options →</td></tr> </table>	HISTORY	History Summary →	Print System History →	Print RO History →	Print Milli-Q History →	Print Options →
MANAGER MENU															
Date and Time →															
Set Points →															
Units →															
Setup →															
User Parameters →															
History →															
Lab closed →															
HISTORY															
History Summary →															
Print System History →															
Print RO History →															
Print Milli-Q History →															
Print Options →															

Item	Description
History Summary	Used to see the day by day history of the Water System.
Print System History	See the section "Printing" for more information.
Print RO History	
Print Milli-Q® History	
Print Options	

Ready Mode

General information

Purpose In READY Mode, water can be dispensed from the POD Unit. The System should be left in READY Mode most of the time.

Display



STANDBY Mode from READY Mode

Display	Action	Result
	Press	

READY Mode – water quality values

The READY Mode screen display is explained below.

READY Mode screen	Explanation
	<p>In this example,</p> <ul style="list-style-type: none"> The permeate water filling the tank has a conductivity of 6 µS/cm. The product water dispensed from the POD Unit has: <ul style="list-style-type: none"> – a resistivity of 18.2 MΩ.cm, – is temperature compensated (TC) at 25°C, and – the TOC value is 4ppb.
	<p>In this example, the System is not dispensing or recirculating water.</p>

Description of Ready Menu

Water Quality

Diagram 1	Diagram 2														
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> <tr><td> </td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →		<table border="1"> <tr><td>WATER QUALITY</td></tr> <tr><td>Permeate Water Quality →</td></tr> <tr><td>Tank Level: 80.0 %</td></tr> <tr><td>Milli-Q Water Quality →</td></tr> <tr><td> </td></tr> </table>	WATER QUALITY	Permeate Water Quality →	Tank Level: 80.0 %	Milli-Q Water Quality →	
READY MENU															
Water Quality →															
Print Menu →															
View Operation →															
Consumables Status →															
Call Millipore →															
Service Tracking →															
Information →															
WATER QUALITY															
Permeate Water Quality →															
Tank Level: 80.0 %															
Milli-Q Water Quality →															

Item	Description
Permeate Water Quality	View the quality of the water filling the Reservoir.
Tank Level	View the level of water in the Reservoir.
MQ Prod Quality	View the quality of water obtained from the POD Unit.

Print Menu

Diagram 1	Diagram 2												
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> <tr><td> </td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →		<table border="1"> <tr><td>PRINT MENU</td></tr> <tr><td>MQ Instant Quality →</td></tr> <tr><td> </td></tr> </table>	PRINT MENU	MQ Instant Quality →	
READY MENU													
Water Quality →													
Print Menu →													
View Operation →													
Consumables Status →													
Call Millipore →													
Service Tracking →													
Information →													
PRINT MENU													
MQ Instant Quality →													

Item	Description
MQ Instant Quality	Print the parameters related to the quality of water delivered from the remote Q-POD® Unit if installed.

Continued on next page

Description of Ready Menu, Continued

View Operation

Diagram 1	Diagram 2															
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> <tr><td> </td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →		<table border="1"> <tr><td>VIEW OPERATION</td></tr> <tr><td>System Operation →</td></tr> <tr><td>System Alerts →</td></tr> <tr><td>System Alarms →</td></tr> <tr><td>System Measures →</td></tr> <tr><td> </td></tr> </table>	VIEW OPERATION	System Operation →	System Alerts →	System Alarms →	System Measures →	
READY MENU																
Water Quality →																
Print Menu →																
View Operation →																
Consumables Status →																
Call Millipore →																
Service Tracking →																
Information →																
VIEW OPERATION																
System Operation →																
System Alerts →																
System Alarms →																
System Measures →																

Item	Description
System Operation	View operating parameters: <ul style="list-style-type: none"> • operating mode, • status of pumps, and • status of UV Lamps.
System Alerts	View a list of active Alert messages. See the Alert Chapter for more information.
System Alarms	View a list of active Alarm messages. See the Alarm Chapter for more information.
System Measures	View: <ul style="list-style-type: none"> • accumulated production time, • pumps electrical data, and • UV Lamps electrical data.

Continued on next page

Description of Ready Menu, Continued

Consumables Status

Diagram 1	Diagram 2																
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →	<table border="1"> <tr><td>CONSUMABLES STATUS</td></tr> <tr><td>Pretreatment →</td></tr> <tr><td>Progard →</td></tr> <tr><td>UV 254 nm Lamp →</td></tr> <tr><td>ASM UV Lamp →</td></tr> <tr><td>UV 185 nm Lamp →</td></tr> <tr><td>Q-Pak →</td></tr> <tr><td>POD Pak →</td></tr> </table>	CONSUMABLES STATUS	Pretreatment →	Progard →	UV 254 nm Lamp →	ASM UV Lamp →	UV 185 nm Lamp →	Q-Pak →	POD Pak →
READY MENU																	
Water Quality →																	
Print Menu →																	
View Operation →																	
Consumables Status →																	
Call Millipore →																	
Service Tracking →																	
Information →																	
CONSUMABLES STATUS																	
Pretreatment →																	
Progard →																	
UV 254 nm Lamp →																	
ASM UV Lamp →																	
UV 185 nm Lamp →																	
Q-Pak →																	
POD Pak →																	

Consumable	Description
Pretreatment	View information about various consumable items. Information may include: <ul style="list-style-type: none"> • installation date, • lifetime remaining, • volume processed, • catalogue number, and • serial number <p>NOTE: Not all of this information is shown for each type of consumable item.</p>
Progard®	
UV 254 nm Lamp	
ASM UV Lamp	
UV 185 nm Lamp	
Q-Pak®	
POD Pak	

Call Millipore SAS

Diagram 1	Diagram 2													
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →	<table border="1"> <tr><td>CALL MILLIPORE</td></tr> <tr><td>Application Specialist →</td></tr> <tr><td>Service Engineer →</td></tr> <tr><td>Tech Service →</td></tr> <tr><td>Other →</td></tr> </table>	CALL MILLIPORE	Application Specialist →	Service Engineer →	Tech Service →	Other →
READY MENU														
Water Quality →														
Print Menu →														
View Operation →														
Consumables Status →														
Call Millipore →														
Service Tracking →														
Information →														
CALL MILLIPORE														
Application Specialist →														
Service Engineer →														
Tech Service →														
Other →														

Item	Description
Application Specialist	View: <ul style="list-style-type: none"> • name, • phone number, and • email address of a Millipore SAS Representative. <p>NOTE: This information is entered by a Millipore SAS Service Representative.</p>
Service Engineer	
Tech Service	
Other	

Continued on next page

Description of Ready Menu, Continued

Service Tracking

Diagram 1	Diagram 2																		
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>InFormation →</td></tr> <tr><td> </td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	InFormation →		<table border="1"> <tr><td>SERVICE TRACKING</td></tr> <tr><td>Installation →</td></tr> <tr><td>Repair →</td></tr> <tr><td>Service Contract →</td></tr> <tr><td>Contract Expires →</td></tr> <tr><td>Next Service →</td></tr> <tr><td>Next Calibration →</td></tr> <tr><td>Next Qualification →</td></tr> <tr><td> </td></tr> </table>	SERVICE TRACKING	Installation →	Repair →	Service Contract →	Contract Expires →	Next Service →	Next Calibration →	Next Qualification →	
READY MENU																			
Water Quality →																			
Print Menu →																			
View Operation →																			
Consumables Status →																			
Call Millipore →																			
Service Tracking →																			
InFormation →																			
SERVICE TRACKING																			
Installation →																			
Repair →																			
Service Contract →																			
Contract Expires →																			
Next Service →																			
Next Calibration →																			
Next Qualification →																			

Item	Description
Installation	View information that was inputted into the System at time of servicing.
Repair	
Service Contract	View information related to upcoming service.
Contract Expires	
Next Service	NOTE: This information is entered by a Millipore SAS Representative.
Next Calibration	
Next Qualification	

Information

Diagram 1	Diagram 2														
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>Print Menu →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>InFormation →</td></tr> <tr><td> </td></tr> </table>	READY MENU	Water Quality →	Print Menu →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	InFormation →		<table border="1"> <tr><td>INFORMATION</td></tr> <tr><td>Flow Schematic →</td></tr> <tr><td>Version →</td></tr> <tr><td>System Information →</td></tr> <tr><td> </td></tr> </table>	INFORMATION	Flow Schematic →	Version →	System Information →	
READY MENU															
Water Quality →															
Print Menu →															
View Operation →															
Consumables Status →															
Call Millipore →															
Service Tracking →															
InFormation →															
INFORMATION															
Flow Schematic →															
Version →															
System Information →															

Item	Description
Flow Schematic	View information that explains the purpose of the major components.
Version	View Software versions.
System Information	View: <ul style="list-style-type: none"> • System Type, • Catalogue Number, • Serial Number, • Installation Date, and • Manufacturing Date.

Continued on next page

Description of Ready Menu, Continued

TOC Curve Check

Diagram 1	Diagram 2														
<table border="1"> <tr><td>READY MENU</td></tr> <tr><td>Water Quality →</td></tr> <tr><td>View Operation →</td></tr> <tr><td>Consumables Status →</td></tr> <tr><td>Call Millipore →</td></tr> <tr><td>Service Tracking →</td></tr> <tr><td>Information →</td></tr> <tr style="background-color: #cccccc;"><td>TOC Curve Check →</td></tr> <tr><td> </td></tr> </table>	READY MENU	Water Quality →	View Operation →	Consumables Status →	Call Millipore →	Service Tracking →	Information →	TOC Curve Check →		<table border="1"> <tr><td>TOC CURVE CHECK</td></tr> <tr><td>Press ✓ to start TOC curve check operation.</td></tr> <tr><td>Press ← to exit.</td></tr> <tr><td> </td></tr> <tr><td> </td></tr> </table>	TOC CURVE CHECK	Press ✓ to start TOC curve check operation.	Press ← to exit.		
READY MENU															
Water Quality →															
View Operation →															
Consumables Status →															
Call Millipore →															
Service Tracking →															
Information →															
TOC Curve Check →															
TOC CURVE CHECK															
Press ✓ to start TOC curve check operation.															
Press ← to exit.															

Item	Description
TOC Curve Check	The TOC Curve Check is used to determine data that is used for the TOC Indicator. A manual TOC Curve Check can be initiated here. Otherwise, the TOC Curve Check is automatically done once per day.

Using the Milli-Q® Direct System

Overview

Introduction

The purpose of this chapter is to explain:

- various ways that water can be dispensed from the System, and
 - how to view information, operating parameters and other things about the System.
-

Contents


This chapter contains the following topics:

Topic	See Page
Dispensing water	64
Viewing water quality	66

Dispensing water

Using the POD Plunger

To dispense water, press down on the POD Unit plunger while in READY Mode.



Position	Water flow
L	Low Flow (push slightly)
M	Medium Flow (push slightly)
H	High Flow (push down and hold, release when done)
H	Continuous high flow (push down and release; push down again to stop).

Continued on next page

Dispensing water, Continued

Volumetric dispensing

Follow the steps below to volumetrically dispense from the POD Unit.

Step	Action	Diagram
1	Make sure the System is in READY Mode.	
2	<ul style="list-style-type: none"> Select Volume. Press 	
3	<ul style="list-style-type: none"> Adjust the volume of water to the required value using and . Press . 	
4	When the volumetric dispensing is finished, the System recirculates water for 3 minutes.	
5	The System stops recirculating water.	

Viewing water quality

Procedure

Follow the steps below to view the water quality.

Step	Action	Diagram
1	Make sure the System is in READY Mode.	
2	<ul style="list-style-type: none"> • Select Menu. • Press 	
3	Select Water Quality. Press .	
4	<ul style="list-style-type: none"> • Select the Water Quality to be viewed. • Press 	
<p>Note The term 'TC' means that the resistivity value is temperature compensated.</p>		
5	Press 3 times on .	

Viewing Operation

Introduction

VIEW OPERATION allows you to see the status of major components. Under the View Operation menu, the following items can be selected:

- System Operation,
- System Alerts,
- System Alarms, and
- System Measures.

System Operation


Follow the steps below to go to the System Operation menu.

Step	Action	Diagram
1	Start in READY Mode.	
2	Select Menu. Press	
3	Select View Operation. Press	
4	Select System Operation. Press	
5	To see more, press	


Continued on next page

Viewing Operation, Continued


System Alerts

<p>An example Alert is shown here. This is an Alert that is currently being displayed on the bottom of the Main Display in READY Mode or in STANDBY Mode.</p>	
<p>When the timer for the UV 185 nm Lamp is reset, then this Alert is no longer shown on the SYSTEM ALERTS LCD.</p>	

System Alarms

<p>An example Alarm is shown here. This is an Alarm that is currently displayed on the Main Display unless you override the display for one hour.</p>	
<p>When the cause of this Alarm is fixed, then this Alarm is no longer shown on the SYSTEM ALARMS LCD.</p>	




System Measures

<p>Various measurements related to the System are shown here.</p>	
---	--

Viewing Consumable Status

Introduction Consumables Status allows you to see information related to the various consumables.

Procedure Follow the steps below to view Consumables Status.

Step	Action	Diagram
1	Start in READY Mode.	
2	Select Menu. Press  .	
3	Select Consumables Status. Press  .	
4	Select the consumable that you would like to see information about. Example The Progard® Cartridge status is shown here. Choose other consumables to see their status.	


Calling Millipore SAS

Introduction

Call Millipore SAS allows you to see contact information.
A Millipore SAS Representative can enter this information into the System.

Procedure

Follow the steps below to view information under Call Millipore SAS.

Step	Action	Diagram
1	Start in READY Mode.	 <p>READY 01 juil. 2009 16:36 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 MΩcm TC TOC: 4 ppb</p>
2	Select Menu. Press  .	 <p>READY MENU Water Quality → Print Menu → View Operation → Consumables Status → Call Millipore → Service Tracking → InFormation →</p>
3	Select Call Millipore SAS. Press  .	 <p>CALL MILLIPORE Application Specialist → Service Engineer → Tech Service → Other →</p>
4	Select the type of Millipore SAS Representative you wish to contact. Press  .	 <p>SERVICE ENGINEER Name: John SMITH Tel: +61 98 9999 Email: John_Smith@Millipore.com ←</p>

Viewing Information

Introduction INFORMATION allows you to view:

- flow schematic information,
- version information, and
- serial number and other information.

Procedure Follow the steps below to see information about the System.

Step	Action	Diagram
1	Start in READY Mode.	
2	Select Menu. Press	
3	Select Information. Press	
4	Select the type of information you wish to view. Two examples are shown below. Press	

Continued on next page

Viewing Information, Continued

Version

The various firmware versions for the System are shown here.

This LCD shows the version used for various components inside the System.	<table border="1"><thead><tr><th data-bbox="938 315 1200 344">VERSION</th></tr></thead><tbody><tr><td data-bbox="938 344 1200 374">Boot Loader: v 1.02</td></tr><tr><td data-bbox="938 374 1200 403">System: v7</td></tr><tr><td data-bbox="938 403 1200 432">EPLD: v1.0</td></tr><tr><td data-bbox="938 432 1200 461">Measure: v1.0</td></tr><tr><td data-bbox="938 461 1200 490">Power Supply: v1.0</td></tr><tr><td data-bbox="938 490 1200 519">POD: v1.0</td></tr><tr><td data-bbox="938 519 1200 548">Tag Reader 1: v1</td></tr></tbody></table>	VERSION	Boot Loader: v 1.02	System: v7	EPLD: v1.0	Measure: v1.0	Power Supply: v1.0	POD: v1.0	Tag Reader 1: v1
VERSION									
Boot Loader: v 1.02									
System: v7									
EPLD: v1.0									
Measure: v1.0									
Power Supply: v1.0									
POD: v1.0									
Tag Reader 1: v1									

System Information

The Catalogue Number, Serial Number and other information are shown here. The Serial Number is something you should reference when you contact Millipore SAS.

This LCD shows information such as the Serial Number and the Catalogue Number. NOTE: The Inst Date (Installation Date) needs to be entered by a Millipore SAS Service Representative. The date is not automatically generated by the System.	<table border="1"><thead><tr><th data-bbox="970 784 1232 813">SYSTEM INFORMATION</th></tr></thead><tbody><tr><td data-bbox="970 813 1232 842">Milli-Q Direct 8</td></tr><tr><td data-bbox="970 842 1232 871">Cat N°: ZR0000000</td></tr><tr><td data-bbox="970 871 1232 900">Serial N°: F6DN27327B</td></tr><tr><td data-bbox="970 900 1232 929">MFG Date: 1 April 2006</td></tr><tr><td data-bbox="970 929 1232 958">Inst Date: 1 June 2006 ←</td></tr></tbody></table>	SYSTEM INFORMATION	Milli-Q Direct 8	Cat N°: ZR0000000	Serial N°: F6DN27327B	MFG Date: 1 April 2006	Inst Date: 1 June 2006 ←
SYSTEM INFORMATION							
Milli-Q Direct 8							
Cat N°: ZR0000000							
Serial N°: F6DN27327B							
MFG Date: 1 April 2006							
Inst Date: 1 June 2006 ←							

Maintenance

Overview

Introduction The purpose of this chapter is to explain the common maintenance needed for a System.

Contents This chapter contains the following topics:

Topic	See Page
Maintenance Schedule	74
Replacing the Progard® Cartridge and Vent Filter	75
Replacing the Q-Pak® Pack	78
Replacing a POD Pak	82
TOC Curve Check	85
Sanitising the RO Cartridge(s)	87
Cleaning the RO Cartridge(s)	90
Cleaning the Inlet Strainer	92
Calibrating the Flow rate	95

Maintenance Schedule

Consumables

Item	Maintenance needed	When
Progard® Cartridge	Replacement	Prompted to by an LCD message.
Q-Pak® Pack	Replacement	
POD Pak	Replacement	Prompted to by an LCD message or as necessary.

UV Lamps

Item	Maintenance needed	When
UV 254 nm Lamp	Replacement	Prompted to by an LCD message.
UV 185 nm Lamp		

NOTE:

It is recommended to have a Millipore SAS Service Representative change the UV Lamps in the system.

The replacement of this lamp involves removing the cover of the system. The instructions for replacing these lamps are not included in this User Manual. The instructions are included with the replacement lamp.

Cleaning/ Sanitisation

Item	Maintenance needed	When
Inlet Strainer	Cleaning	Prompted to by an LCD message or as necessary.
RO Cartridge(s)	Cl ₂ cleaning	When prompted to by an LCD message.
	pH Cleaning	As necessary.
System	Sanitisation	Contact Millipore SAS for more details.

Calibrating the flow rate

Item	Maintenance needed	When
Flow rate	Recalibration	New Consumable, Sensor or change to Feed water. See 'Calibrating the flow rate' for more information.

TOC Curve Check

Item	Maintenance needed	When
TOC Indicator	Update TOC Curve Check	New Q-Pak® Pack installed, or when prompted to by an LCD message.

Replacing the Progard® Cartridge and Vent Filter

When

The Progard® Cartridge and Tank Vent Filter should be replaced when the following Alert message is displayed.

- Alert message = REPLACE PROGARD CARTRIDGE AND TANK VENT FILTER
-



Attention


The Progard® Cartridge must be flushed after it is installed.

Continued on next page

Replacing the Progard® Cartridge and Vent Filter, Continued

Removing

Follow the steps below to remove the used Progard® Cartridge.


Step	Action	Diagram
1	Place the System into STANDBY Mode.	
2	Push the POD Plunger down once to depressurise the System. After water stops being dispensed, push down the POD Plunger again.	
3	Open the System right door. Remove the used Progard® Cartridge.	
4	In a few moments, the System indicates that the Progard® Cartridge is removed.	

Continued on next page

Replacing the Progard® Cartridge and Vent Filter, Continued

Placing

Follow the steps below to install a new Progard® Cartridge.



Step	Action	Diagram								
1	Remove the covers on the 2 ports of the Progard® Cartridge. Wet the O-rings with water.									
2	Install the Progard® Cartridge until it is fully seated. Close the right door.									
3	When a new Progard® Cartridge is installed, the LCD looks like this.	<table border="1"> <tr> <td>INSTALL PROGARD</td> </tr> <tr> <td>A new Progard has been installed.</td> </tr> <tr> <td>Catalogue N°: PR0G000T3</td> </tr> <tr> <td>Lot N°: F8DN27324.</td> </tr> <tr> <td>Press → to start Progard Flush/cleaning.</td> </tr> </table>	INSTALL PROGARD	A new Progard has been installed.	Catalogue N°: PR0G000T3	Lot N°: F8DN27324.	Press → to start Progard Flush/cleaning.			
INSTALL PROGARD										
A new Progard has been installed.										
Catalogue N°: PR0G000T3										
Lot N°: F8DN27324.										
Press → to start Progard Flush/cleaning.										
4	Press  .	<table border="1"> <tr> <td>INSTALL PROGARD</td> </tr> <tr> <td>Progard Flush procedure in progress.</td> </tr> <tr> <td>Remaining Time: XX min.</td> </tr> <tr> <td>Press → to cancel.</td> </tr> </table>	INSTALL PROGARD	Progard Flush procedure in progress.	Remaining Time: XX min.	Press → to cancel.				
INSTALL PROGARD										
Progard Flush procedure in progress.										
Remaining Time: XX min.										
Press → to cancel.										
5	When the Progard® Cartridge flush has finished, the Water System goes to READY Mode.	<table border="1"> <tr> <td>READY</td> </tr> <tr> <td>02 juil. 2009 11:45</td> </tr> <tr> <td>Menu →</td> </tr> <tr> <td>Tank: Standby →</td> </tr> <tr> <td>0 % Volume →</td> </tr> <tr> <td>Perm C: 6.0 µS/cm TC</td> </tr> <tr> <td>MQ Res: - - - MΩcm TC</td> </tr> <tr> <td>TOC: - - - ppb</td> </tr> </table>	READY	02 juil. 2009 11:45	Menu →	Tank: Standby →	0 % Volume →	Perm C: 6.0 µS/cm TC	MQ Res: - - - MΩcm TC	TOC: - - - ppb
READY										
02 juil. 2009 11:45										
Menu →										
Tank: Standby →										
0 % Volume →										
Perm C: 6.0 µS/cm TC										
MQ Res: - - - MΩcm TC										
TOC: - - - ppb										

Replacing the Q-Pak® Pack

When The Q-Pak® Pack should be replaced when one of the following Alarm or Alert messages is displayed.

- Alarm message = MILLI-Q RES < SP, REPLACE Q-PAK
- Alert message = REPLACE Q-PAK PACK

Removing Remove the used Q-Pak® Pack by following the steps below.

Step	Action	Diagram
1	Place the system into STANDBY Mode.	
2	Push the POD Plunger down once to depressurise the System. After water stops being dispensed, push down the POD Plunger again.	
3	Open the System left door. Lift up the Pack Locking Handle.	

Continued on next page


Replacing the Q-Pak® Pack, Continued

Removing (continued)

Step	Action	Diagram
4	Remove the used Q-Pak® Pack.	
5	The System will indicate that the Q-Pak® Pack is removed in a few moments.	

Placing

Follow the steps below to install a new Q-Pak® Pack.

Step	Action	Diagram
1	Remove the covers on the 2 ports of the Q-Pak® Pack. Look inside the ports. Make sure the rubber O-rings are firmly in place. Wet the O-rings with water.	
2	Push the top of the Q-Pak® Pack into the ports on the System. Push on the bottom of the Q-Pak® Pack.	

Continued on next page

Replacing the Q-Pak® Pack, Continued

Placing (continued)

Step	Action	Diagram	
3	Push the Pack Locking Handle down. Close the left door.		
4	When a new Q-Pak® Pack is installed, the screen looks like this. Press  .	<table border="1"> <tr> <td> INSTALL Q-PAK A new Q-PAK has been installed. Catalogue N°: QPAK00TEX Lot N°: F6DN27325. ← </td> </tr> </table>	INSTALL Q-PAK A new Q-PAK has been installed. Catalogue N°: QPAK00TEX Lot N°: F6DN27325. ←
INSTALL Q-PAK A new Q-PAK has been installed. Catalogue N°: QPAK00TEX Lot N°: F6DN27325. ←			

Rinsing


The Q-Pak® pack, when newly installed, needs to be rinsed. This ensures optimal water quality.

Step	Action	Diagram
1	Locate the clear tubing and the barbed fitting from the System accessories bag. Screw the barbed fitting onto the POD Unit. NOTE: Do not use any white tape on the threads of the barbed fitting. An O-ring is located inside the POD Unit. Push one end of the clear tubing onto the end of the barbed fitting. Place the other end of the clear tubing into a sink.	

Continued on next page

Replacing the Q-Pak® Pack, Continued

Rinsing (continued)

Step	Action	Diagram
2	The System must be in READY Mode.	
3	Push the plunger down on the POD Unit.	
4	Dispense water for about 10 minutes. This flushes out any trapped air in most of the System. This also rinses off the purification media located in the Q-Pak® Pack.	
5	Leave the System in READY Mode when finished.	

Manual TOC Curve Check

A TOC Curve Check should be performed when the Q-Pak® Pack has been replaced. Refer to the TOC Curve Check section for more information.

Flow Rate Calibration

The volumetric dispensing flow rate should be calibrated when the Q-Pak® Pack has been replaced. Refer to the procedure in the Installation chapter.

Replacing a POD Pak

Basing on flow rate

One possible reason for a decrease in Milli-Q® Water flow rate is a clogged POD Pak. The POD Pak should be replaced when it appears to be clogged. For Millipak® and BioPak® final filters, make sure the POD Pak is not air-locked. Dispense water and open the vent to see if there is any trapped air. Close the vent after this.

Basing on LCD message

The POD Pak needs replacement when the following Alert message is displayed.


- Alert message = REPLACE POD PAK

Placing and flushing

Follow the instructions delivered with the POD Pak.

Registering

The POD Pak installation has to be registered. Follow the steps below to register the installation of the POD Pak.

Step	Action	Diagram
1	Start in STANDBY Mode.	
2	Select Menu. Press  .	
3	Select Maintenance. Press  .	

Continued on next page

Replacing a POD Pak, Continued

Registering (continued)

Step	Action	Diagram
4	Scroll down to Install POD Pak by pressing  .	 <p>MAINTENANCE Clean Strainer → Install Progard → Install new RO → Install UV 254 nm Lamp → Install UV 185 nm Lamp → Install Q-Pak → Install POD Pak →</p>
5	Press  .	 <p>INSTALL POD PAK</p>
6	Press  .	 <p>INSTALL POD PAK Select the POD Pak that you wish to install. Press → to continue or ← to exit.</p>
7	In this example, the replacement POD Pak is a Millipak®. Press  .	 <p>INSTALL POD PAK Millipak → BioPak → EDS-Pak → Other Pod Pak A → Other Pod Pak B → No Filter →</p>
8	Press  .	 <p>INSTALL POD PAK Follow the instructions delivered with the new POD Pak and press √. ←</p>

Continued on next page

Replacing a POD Pak, Continued

Registering (continued)

Step	Action	Diagram
9	Press  .	
10	Press 3 times on  .	

Flow Rate Calibration

The volumetric dispensing flow rate should be calibrated when a POD Pak has been replaced. Refer to the procedure in the Installation chapter.

TOC Curve Check

When A TOC Curve Check should be done when the Q-Pak® pack is replaced or when TOC values are fluctuating to ensure that the TOC values displayed are valid.

NOTE:

The TOC Curve Check can be started manually following the replacement and flushing of the Q-Pak® pack. The System automatically performs a TOC Curve Check once per day. So, if the TOC Curve Check is not started manually following the consumables replacement, it will be done automatically within the next 24 hours or sooner.

Procedure Follow the steps below to perform a TOC Curve Check.

Step	Action	Diagram
1	Go to READY Mode.	
2	Select Menu. Press	
3	Select TOC Curve Check. Press	

Continued on next page

TOC Curve Check, Continued

Procedure (continued)

Step	Action	Diagram
4	Press  .	<div data-bbox="970 342 1230 607" style="border: 1px solid black; padding: 5px;"> <p>TOC CURVE CHECK</p> <p>The system is now in TOC curve check processing.</p> <p>Task Completion: XX min</p> <p>Press ← to cancel and exit.</p> </div>
5	After approximately 10 minutes, the System returns to READY Mode.	<div data-bbox="970 622 1230 887" style="border: 1px solid black; padding: 5px;"> <p>READY</p> <p>16 juil. 2009 16:43</p> <p style="text-align: right;">Menu →</p> <p>Tank: Standby →</p> <p>80.0 % Volume →</p> <p>Perm C: 6.0 µS/cm TC</p> <p>MQ Res: 18.2 MΩcm TC</p> <p>TOC: 4 ppb</p> </div>

Sanitising the RO Cartridge(s)

When

RO cartridge sanitisation is required to prevent bacteria development in the cartridge. To maintain optimum RO performance, perform sanitisation when the following alert message is displayed:

- Alert message = PERFORM RO CL2 CLEANING.



Wear Eye Safety Glasses and Laboratory Gloves and other appropriate safety equipment when sanitising the RO Cartridge(s).

Opening the Sanitisation Port

Follow the steps below to open the Sanitisation Port.

Step	Action	Diagram
1	<ul style="list-style-type: none"> • Go to STANDBY Mode. • Allow the Water System to depressurise for a few seconds. 	
2	Use the Sanitisation Port Removal Tool and loosen the cap.	
3	Remove the cap.	

Closing the Sanitisation Port

Reverse the steps above.

NOTE:




Do not use the Sanitisation Port Removal Tool to tighten the cap.

Continued on next page

Sanitising the RO Cartridge(s), Continued

Sanitising



Follow the steps below to sanitise the RO Cartridge(s).

Step	Action	Diagram
1	Place a chlorine tablet into the Sanitisation Port.	
2	Put the cap back on and hand-tighten it.	
3	Go to STANDBY Mode.	
4	<ul style="list-style-type: none"> • Select Menu. • Press . 	

Continued on next page

Sanitising the RO Cartridge(s), Continued

Sanitising (continued)

Step	Action	Diagram
5	<ul style="list-style-type: none"> • Select Sanitise/Clean. • Press . 	
6	<ul style="list-style-type: none"> • Select RO CL2 Cleaning. • Press . 	
7	<ul style="list-style-type: none"> • Press . • The RO CL2 cleaning mode will last 19 minutes. 	
8	When the cleaning is finished, the Water System automatically goes into READY Mode.	

Continued on next page



Cleaning the RO Cartridge(s)

When Cleaning the RO cartridge may be required when feed water contains excess of some chemicals.
Perform cleaning if required, after recommendation from your Millipore SAS Service Representative.



Wear Eye Safety Glasses and Laboratory Gloves and other appropriate safety equipment when cleaning the RO Cartridge(s).




Cleaning Follow the steps below to clean the RO Cartridge(s).

Step	Action	Diagram
1	Open the sanitisation port as described in the previous section.	
2	Place a cleaning agent pouch (ROClean™ A or ROClean™ B) into the Sanitisation Port.	
3	Put the cap back on and hand-tighten it. NOTE: The chemical in the pouch will dissolve during the pH Cleaning sequence. Remove the empty pouch the next time the Sanitisation Port cap is removed.	
4	Go to STANDBY Mode.	

Continued on next page

Cleaning the RO Cartridge(s), Continued

Cleaning (continued)

Step	Action	Diagram
5	<ul style="list-style-type: none"> • Select Menu. • Press . 	
6	<ul style="list-style-type: none"> • Select Sanitise/Clean. • Press . 	
7	<ul style="list-style-type: none"> • Select RO pH Cleaning. • Press . 	
8	<ul style="list-style-type: none"> • Press . • The RO pH cleaning will last 142 minutes. 	
9	When the pH Cleaning is finished, the Water System automatically returns to READY Mode.	

Cleaning the Inlet Strainer

- Purpose**
- The purpose of the Inlet Strainer is to prevent a large particle from entering the System.
 - If the Inlet Strainer becomes clogged, then feed water does not flow freely to the System.
 - Cleaning the Inlet Strainer removes any trapped debris.
-

When

The Inlet Strainer should be cleaned when the following Alert message is displayed.

- Alert message = EXAMINE INLET STRAINER

The Inlet Strainer should also be cleaned whenever you suspect it is clogged.

Procedure

Follow the steps below to clean the Inlet Strainer.




Step	Action
1	Go to STANDBY Mode.
2	Shut off the feed water supply.
3	Unscrew the Inlet Strainer from the feed water supply.
4	Detach the tubing on the other end of the Inlet Strainer.
5	Flush water backwards through the Inlet Strainer.
6	Apply 3 to 4 turns of new white tape to the threads of the feed water pipe.
7	Screw the Inlet Strainer back onto the feed water pipe.
8	Attach the tubing to the other end of the Inlet Strainer.
9	Open the feed water supply valve.
10	Go to READY Mode.

Continued on next page

Cleaning the Inlet Strainer, Continued

Registering

Follow the steps below to register the cleaning of the Inlet Strainer.

Step	Action	Diagram
1	Go to STANDBY Mode.	 <p>STANDBY 15 Jul. 2009 16:23 Menu → Ready →</p>
2	Select Menu. Press  .	 <p>STANDBY MENU Maintenance → Sanitise/Clean → Language → Manager Menu →</p>
3	Select Maintenance. Press  .	 <p>MAINTENANCE Clean Strainer → Install Progard → Install new RO → Install UV 254 nm Lamp → Install UV 185 nm Lamp → Install Q-Pak → Install POD Pak →</p>
4	Select Clean Strainer. Press  .	 <p>CLEAN STRAINER</p>
5	A picture is shown. Press  .	 <p>CLEAN STRAINER See Maintenance Chapter in the User Manual For more information. Press ✓ after cleaning or ← to exit.</p>

Continued on next page

Cleaning the Inlet Strainer, Continued

Procedure
(continued)

Step	Action	Diagram
6	Press  .	<div data-bbox="970 353 1230 616" style="border: 1px solid black; padding: 5px;"> <p>CLEAN STRAINER</p> <p>The strainer cleaning date is registered. Next maintenance in 365 days. Press ← to exit.</p> </div>
7	Press 3 times on  .	<div data-bbox="970 645 1230 907" style="border: 1px solid black; padding: 5px;"> <p>STANDBY</p> <p>01 juil. 2009 16:55</p> <p style="text-align: right;">Menu →</p> <p style="text-align: right;">Ready →</p> </div>
8	Go to READY Mode.	<div data-bbox="970 925 1230 1187" style="border: 1px solid black; padding: 5px;"> <p>READY</p> <p>01 juil. 2009 16:55</p> <p style="text-align: right;">Menu →</p> <p>Tank: Standby →</p> <p>80.0 % Volume →</p> <p>Perm C: 6.0 µS/cm TC</p> <p>MQ Res: 18.2 MΩ·cm TC</p> <p>TOC: 4 ppb</p> </div>

Calibrating the Flow rate

When

The flow rate should be calibrated when:

- a new consumable is installed such as:
 - POD Pak, or
 - Q-Pak® Pack,
 - a sensor or major component is changed,
 - volumetric dispensing is not accurate.
-

Procedure

Follow the procedure shown in the Installation Chapter.

Alarms

Overview

Introduction

The purpose of this chapter is to explain the Alarm messages shown on a System. Specifically, this chapter explains how:

- an Alarm message is displayed,
 - to read an Alarm message,
 - to cancel an Alarm, and
 - a list of Alarm messages is shown.
-

Contents

This chapter contains the following topics:

Topic	See Page
Alarm Information	97
Summary of Alarm messages	101

Alarm Information

Definition An Alarm message is a way of informing you that immediate attention is required for the System.



Alarm shown – what to do?

It is not recommended to use the System when an Alarm message is shown. Contact Millipore SAS if an Alarm message is shown and the problem can not be resolved.

Types The following table summarizes the different types of Alarm messages.

Type	Description
Alarm stop	Some Alarms automatically stop the System from dispensing water. An example of this is the Alarm message PROGARD CARTRIDGE OUT. The text display of this type of Alarm can be cancelled for one hour by using the Keypad.
Alarm	Some Alarms do not automatically stop the System from dispensing water. An example of this is the Alarm message MILLI-Q T < MIN. The text display of this type of Alarm can be cancelled for one hour by using the Keypad.

Main Display The Alarm message is shown superimposed on the Main Display. The red LED is lit steadily when an Alarm message is shown. In this example, the Alarm Message MILLI-Q T > MAX is shown.




Continued on next page

Alarm Information, Continued

System Alarms

When an Alarm is shown, it is listed under the System Alarms LCD. See the section <View Operation> for information on how to access this LCD.



Viewing an Alarm Message

Follow the steps below to view an Alarm message.

Step	Action	Diagram
1	The Alarm message is shown superimposed on the Main Display.	<p>The diagram shows the main display with the alarm message 'MILLI-Q T > MAX' superimposed. The main display text includes 'READY', '02 Jul, 2009 15:07', 'Tank', '80.0', 'Perm', 'MQ RPRESS →', and 'TOC: 4 ppb'. The alarm message is shown in a larger font and is partially obscured by a black box.</p>
2	Press  .	<p>See Alarms Chapter in the User Manual For more information. Press ✓ to cancel the display of this alarm For one hour or press ← to exit.</p>
3	Press  .	<p>The diagram shows the main display with the alarm message 'MILLI-Q T > MAX' superimposed. The main display text includes 'READY', '02 Jul, 2009 15:07', 'Tank', '80.0', 'Perm', 'MQ RPRESS →', and 'TOC: 4 ppb'. The alarm message is shown in a larger font and is partially obscured by a black box.</p>

Continued on next page

Alarm Information, Continued


Canceling an Alarm message

The display of an Alarm message can be cancelled by:

- fixing the cause of the Alarm, or
- using the Keypad. This cancels the display of the Alarm message for 1 hour.

Alarm – before cancelling

In this example, the Alarm message is MILLI-Q T > MAX.



Canceling an Alarm message procedure

Follow the steps below to cancel an Alarm message.

Step	Action	Diagram
1	The Alarm message is shown superimposed on the Main Display.	<p>READY 02 Jul 2009 15:07 MILLI-Q T > MAX Tank nu → 80.0 by → Perm ve → MQ PRESS → TOC: 4 ppb</p>
2	Press .	<p>See Alarms Chapter in the User Manual For more information. Press ✓ to cancel the display of this alarm For one hour or press ← to exit.</p>
3	Press .	The display of the Alarm is cancelled for one hour. It appears after one hour unless the cause of the Alarm is fixed.

Continued on next page

Alarm Information, Continued

Alarm – after cancelling the text display

Main Display	LEDs	Main Display
<pre> READY 02 juil. 2009 15:08 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 MΩ.cm TC TOC: 4 ppb </pre>		<pre> SYSTEM ALARMS Milli-Q T > Max </pre>

Alarm – fixed

Now suppose a Millipore SAS Service Representative fixes the cause of the Alarm.

Main Display	LEDs	Main Display
<pre> READY 02 juil. 2009 15:08 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 MΩ.cm TC TOC: 4 ppb </pre>		<pre> SYSTEM ALARMS No Alarms </pre>

Summary of Alarm messages

Alarm stop messages

LCD message	What it means
FLOW AUTO STOP	<ul style="list-style-type: none"> • This alarm stops the Milli-Q® portion of the system. • A safety feature of the Milli-Q® Direct system has automatically stopped dispensing water after 40 minutes to avoid risks of overflow. • Push the POD Unit Plunger all the way down and release. • This resets the dispenser timer and makes the POD Unit available for dispensing.
INCORRECT PROGARD CARTRIDGE	<ul style="list-style-type: none"> • This alarm stops the RO portion of the system. • The System does not recognize the type of Progard® Cartridge being installed. • Contact Millipore SAS.
INCORRECT Q-PAK PACK	<ul style="list-style-type: none"> • This alarm stops the Milli-Q® portion of the system. • The System does not recognize the type of Q-Pak® Pack being installed. • Contact Millipore SAS.
LOW FEED WATER PRESSURE	<ul style="list-style-type: none"> • This alarm stops the RO portion of the system. • Check Feed water pressure and rectify. • Go to STANDBY Mode and go to READY Mode to release any trapped air in the Water System. • Contact Millipore SAS if the problem persists.
PERMEATE C > SP	<ul style="list-style-type: none"> • This alarm stops the RO portion of the system. • The Permeate conductivity is above the set point. • Contact Millipore SAS.
POD LOCKED	<ul style="list-style-type: none"> • This alarm stops the Milli-Q® portion of the system. • The POD Unit was left in the open position. • Push the Plunger all the way down and release.

Continued on next page

Summary of Alarm messages, Suite

Alarm stop messages (suite)

LCD message	What it means
PROGARD CARTRIDGE OUT	<ul style="list-style-type: none"> • This alarm stops the RO portion of the system. • The Progard® Cartridge is not installed correctly or it has been removed. • Verify that the Progard® Cartridge is installed correctly. • Contact Millipore SAS if the problem continues.
Q-PAK PACK OUT	<ul style="list-style-type: none"> • This alarm stops the Milli-Q® portion of the system. • The Q-Pak® Pack is not installed correctly or it has been removed. • Verify that the Q-Pak® Pack is installed correctly. • Contact Millipore SAS if the problem continues.
TANK EMPTY	<ul style="list-style-type: none"> • This alarm stops the Milli-Q® portion of the system. • The System has detected an empty Reservoir. • Refill the Reservoir. • Verify that the Reservoir level sensor is plugged into the System Cabinet.
WATER DETECTED	<ul style="list-style-type: none"> • This alarm stops the whole system to prevent risks of flood. • A Water Sensor (an accessory connected to the System) has detected water on the surface where it is located. This may be caused by a leak. • Clean up the spilled water. • Make sure the source of the leak is fixed. • Place the system in Standby Mode, then Ready mode.

Ce sujet continue page suivante

Summary of Alarm messages, Continued

Alarm messages

CLEANING CANCELLED	<ul style="list-style-type: none"> • A cleaning mode was cancelled and was not fully completed. • Go to STANDBY Mode and then go to READY Mode. • The Milli-Q® System will go into a 15 minute FLUSH Mode. The system will then automatically fill the Reservoir.
MILLI-Q INTER R > MAX	<ul style="list-style-type: none"> • The Intermediate resistivity is out of measurement range. • Contact Millipore SAS.
MILLI-Q INTER R < MIN	
MILLI-Q INTER T < MIN	<ul style="list-style-type: none"> • The Intermediate temperature is out of measurement range. • Contact Millipore SAS.
MILLI-Q INTER T > MAX	
MILLI-Q RES < SP, REPLACE Q-PAK	<ul style="list-style-type: none"> • The Milli-Q® Water resistivity is below the set point. • Dispense water to eliminate any trapped air in the System. • If the issue persists, replace the Q-Pak® Pack.
MILLI-Q RES > MAX	<ul style="list-style-type: none"> • The Milli-Q® Water resistivity is out of measurement range. • Contact Millipore SAS.
MILLI-Q T < MIN	<ul style="list-style-type: none"> • The Milli-Q® Water temperature is out of measurement range. • Contact Millipore SAS.
MILLI-Q T > MAX	
MILLI-Q TOC > SP	<ul style="list-style-type: none"> • The TOC is above the set point. • Contact Millipore SAS.
PERMEATE C < MIN	<ul style="list-style-type: none"> • The Permeate conductivity is out of measurement range. • Contact Millipore SAS.
PERMEATE C > MAX	
RO FEED C < MIN	<ul style="list-style-type: none"> • The Feed water conductivity is out of measurement range. • Contact Millipore SAS.
RO FEED C > MAX	
RO FEED T < MIN	<ul style="list-style-type: none"> • The Feed water temperature is out of measurement range. • Contact Millipore SAS.
RO FEED T > MAX	

Alerts

Overview

Introduction

The purpose of this chapter is to explain the Alert messages shown on a System. Specifically, this chapter explains how:

- an Alert message is displayed,
 - to read an Alert message,
 - to cancel an Alert, and
 - a list of Alert messages is shown.
-

Contents

This chapter contains the following topics:

Topic	See Page
Alert information	105
Summary of Alert messages	110

Alert information

Purpose An Alert message corresponds to a maintenance request. Most of the Alert messages are related to the replacement of a consumable.

Types The following table summarizes the different types of Alert messages.

Type	Description
Minor Alert	A minor alert message indicates that a maintenance action is needed within a number of days.
Major Alert	A major Alert message corresponds to an immediate maintenance request.

Examples An example of a minor alert message would be REPLACE POD PAK IN 15 DAYS. An example of a major alert message would be REPLACE POD PAK.

Main Display An Alert message is shown on the bottom of the Main Display. In this example, the Alert message REPLACE POD PAK scrolls across the bottom of the LCD.

READY
08 juil. 2009 10:42
Menu →
Tank: Standby →
80.0 % Volume →
Perm C: 6.0 µS/cm TC
MQ Res: 18.2 MΩ.cm TC
TOC: 4 ppb
* REPLACE POD PAK *** PF

The yellow LED is lit steadily when an Alert message is shown. However, if an Alert and an Alarm are both present, then only the red LED is lit. When an Alert is shown, it is listed under the System Alerts LCD. To access the System Alerts LCD, see the Section View Operation.

SYSTEM ALERTS
Replace POD Pak

Continued on next page

Alert information, Continued

Viewing an Alert Message

Follow the steps below to view an Alert message.

Step	Action	Diagram
1	Alert messages appear on the bottom line of the screen when the system is in READY mode or in STANDBY mode.	
2	Press as many times as required to reach the bottom line.	
3	Press to access detailed information about the alert.	<p>The POD Pak installed on Point of Distribution should be replaced. Please make sure to replace it on time for optimal system performance. See Alerts Chapter in the User Manual for more information.</p>
4	Press to view full information about the alert if needed.	<p>make sure to replace it on time for optimal system performance. See Alerts Chapter in the User Manual for more information. Press to cancel the text display of this alert or press to exit.</p>
5	Press .	

Continued on next page


Alert information, Continued

Canceling a Minor Alert message - procedure

A Minor alert message can be cancelled by:

- performing the maintenance action (i.e. replace consumable),
- using the Keypad (see below), or
- a Major Alert message is shown. This eliminates the Minor Alert message.

Example: Before cancelling, the Minor Alert message is REPLACE POD PAK IN 15 DAYS.



Follow the steps below to cancel a Minor Alert message.

Step	Action	Diagram
1	Press  .	<p>READY 08 Jul. 2009 10:44 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 MΩcm TC TOC: 4 ppb * REPLACE POD PAK IN 15</p>
2	Press  .	<p>The POD Pak installed on Point of Distribution should be replaced in 15 days. Please make sure to replace it on time for optimal system performance. See Alerts Chapter in the User Manual</p>
3	Press  .	The display of the Minor Alert is cancelled.

Continued on next page

Alert information, Continued

Minor Alert – after cancelling

The Alert message has been cancelled but the cause of the message is still active.

Main Display	LEDs	Main Display
<pre> READY 08 juil. 2009 10:44 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 Maxcm TC TOC: 4 ppb </pre>		<pre> SYSTEM ALERTS Replace POD Pak in 15 days </pre>

Minor Alert – consumable replaced

The Alert message has been cancelled when the POD Pak has been replaced.

Main Display	LEDs	Main Display
<pre> READY 08 juil. 2009 10:44 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 Maxcm TC TOC: 4 ppb </pre>		<pre> SYSTEM ALERTS No Alerts </pre>

Cancelling a Major Alert message – procedure

A Major Alert message can be cancelled by:

- performing the maintenance action (i.e. replace consumable), or
- using the Keypad. This cancels the display of the Major Alert message for 24 hours.

Example: Before cancelling, the Major Alert message is REPLACE POD PAK.

Main Display	LEDs	Main Display
<pre> READY 08 juil. 2009 10:45 Menu → Tank: Standby → 80.0 % Volume → Perm C: 6.0 µS/cm TC MQ Res: 18.2 Maxcm TC TOC: 4 ppb * REPLACE POD PAK *** Pf </pre>		<pre> SYSTEM ALERTS Replace POD Pak </pre>

A Major Alert message can be cancelled using the Keypad. This is done in the same way that a Minor Alert message is cancelled.

The display of the Major Alert is cancelled for 24 hours. It appears again after 24 hours unless the maintenance action is performed.

Continued on next page

Alert information, Continued

**Major Alert –
after cancelling**

The Alert message has been cancelled but the cause of the message is still active.

Main Display	LEDs	Main Display
<p>READY</p> <p>08 Jul. 2009 10:45</p> <p>Menu →</p> <p>Tank: Standby →</p> <p>80.0 % Volume →</p> <p>Perm C: 6.0 µS/cm TC</p> <p>MQ Rest: 18.2 MΩcm TC</p> <p>TOC: 4 ppb</p>		<p>SYSTEM ALERTS</p> <p>Replace POD Pak</p>

**Major Alert –
consumable
replaced**

The Alert message has been cancelled when the POD Pak has been replaced.

Main Display	LEDs	Main Display
<p>READY</p> <p>08 Jul. 2009 10:45</p> <p>Menu →</p> <p>Tank: Standby →</p> <p>80.0 % Volume →</p> <p>Perm C: 6.0 µS/cm TC</p> <p>MQ Rest: 18.2 MΩcm TC</p> <p>TOC: 4 ppb</p>		<p>SYSTEM ALERTS</p> <p>No Alerts</p>

Summary of Alert messages

Alert messages

LCD message	What it means
CALIBRATION VISIT OVERDUE XX DAYS	<ul style="list-style-type: none"> The System has determined that a Calibration Visit is overdue. Contact Millipore SAS.
CHECK ASM UV LAMP	<ul style="list-style-type: none"> The ASM UV Lamp is not turning on. Contact Millipore SAS.
CHECK UV 185 NM LAMP	<ul style="list-style-type: none"> The UV 185 nm Lamp is not turning on. Contact Millipore SAS.
CHECK UV 254 NM LAMP	<ul style="list-style-type: none"> The UV 254 nm Lamp is not turning on. Contact Millipore SAS.
EXAMINE INLET STRAINER	<ul style="list-style-type: none"> The System has determined that it is time to clean the Inlet Strainer. Clean the Inlet Strainer and reset the message.
LOW RO PUMP PRESSURE	<ul style="list-style-type: none"> The Water System has determined that the RO Pump Pressure is below specification. Contact Millipore SAS.
NEXT CALIBRATION VISIT IN XX DAYS	<ul style="list-style-type: none"> The System is prompting you that a Calibration Visit should be scheduled. Contact Millipore SAS.
NEXT QUALIFICATION VISIT IN XX DAYS	<ul style="list-style-type: none"> The System is prompting you that a Qualification Visit should be scheduled. Contact Millipore SAS.
NEXT SERVICE VISIT IN XX DAYS	<ul style="list-style-type: none"> The System is prompting you that a Service Visit should be scheduled. Contact Millipore SAS.
NO RESPONSE FROM DHCP SERVER	<ul style="list-style-type: none"> Contact your network administrator. Restart the System.
PERFORM RO CL2 CLEANING	<ul style="list-style-type: none"> The Water System has determined that it is time to perform a RO Cl2 cleaning. Start a RO Cl2 cleaning. The timer is reset automatically.

Continued on next page

Summary of Alert messages, Continued

Alert messages (continued)

LCD message	What it means
PERFORM TOC CURVE CHECK	<ul style="list-style-type: none"> • The System has determined that a TOC Curve Check should be performed. • Perform a TOC Curve Check or wait until one is done automatically. • If this message persists after 24 hours of its appearance, then contact Millipore SAS.
QUALIFICATION VISIT OVERDUE XX DAYS	<ul style="list-style-type: none"> • The System has determined that a Qualification Visit is overdue. • Contact Millipore SAS.
REPLACE ASM UV LAMP	<ul style="list-style-type: none"> • The Water System has determined that the ASM UV Lamp should be replaced. • Contact Millipore SAS.
REPLACE ASM UV LAMP IN XX DAYS	<ul style="list-style-type: none"> • The Water System has determined that the ASM UV Lamp on the Reservoir should be replaced in XX days, where XX is 15, ..., 1. • Contact Millipore SAS.
REPLACE EXTERNAL PRETREATMENT	<ul style="list-style-type: none"> • The Water System has determined that the external pretreatment (optional) should be replaced. • Consult the documentation supplied with the external pretreatment for more information.
REPLACE EXTERNAL PRETREATMENT IN XX DAYS	<ul style="list-style-type: none"> • The Water System has determined that the external pretreatment (optional) should be replaced in XX days, where XX is 15, ..., 1. • Consult the documentation supplied with the external pretreatment for more information.
REPLACE POD PAK	<ul style="list-style-type: none"> • The System has determined that the POD PAK needs replacement. • Replace the POD Pak and reset the timer.
REPLACE POD PAK IN XX DAYS	<ul style="list-style-type: none"> • The System has determined that the POD PAK should be replaced in XX days, where XX is 15, ..., 1. • Replace the POD Pak and reset the timer.

Continued on next page

Summary of Alert messages, Continued

Alert messages (continued)

REPLACE Q-PAK PACK	<ul style="list-style-type: none"> • The System has determined that the Q-Pak® Pack should be replaced. • Replace the Q-Pak® Pack.
REPLACE Q-PAK PACK IN XX DAYS	<ul style="list-style-type: none"> • The System has determined that the Q-Pak® Pack should be replaced in XX days, where XX is 15, ..., 1. • Replace the Q-Pak® Pack.
REPLACE PROGARD AND TANK VENT FILTER	<ul style="list-style-type: none"> • The System has determined that the Progard® Cartridge and the Vent Filter should be replaced. • Replace the Progard® Cartridge and the Vent Filter.
REPLACE PROGARD AND TANK VENT FILTER IN XX DAYS	<ul style="list-style-type: none"> • The System has determined that the Progard® Cartridge and the Vent Filter should be replaced in XX days, where XX is 15, ..., 1. • Replace the Progard® Cartridge and the Vent Filter.
REPLACE UV 185 NM LAMP	<ul style="list-style-type: none"> • The System has determined that the UV 185 nm Lamp should be replaced. • Contact Millipore SAS.
REPLACE UV 185 NM LAMP IN XX DAYS	<ul style="list-style-type: none"> • The System has determined that the UV 185 nm Lamp should be replaced in XX days, where XX is 15, ..., 1. • Contact Millipore SAS.
REPLACE UV 254 NM LAMP	<ul style="list-style-type: none"> • The Water System has determined that the UV 254 nm Lamp should be replaced. • Contact Millipore SAS.
REPLACE UV 254 NM LAMP IN XX DAYS	<ul style="list-style-type: none"> • The Water System has determined that the UV 254 nm Lamp should be replaced in XX days, where XX is 15, ..., 1. • Contact Millipore SAS.
RO REJECTION < SP	<ul style="list-style-type: none"> • The RO % Rejection is below the set point. • Contact Millipore SAS.
SERVICE VISIT OVERDUE XX DAYS	<ul style="list-style-type: none"> • The System has determined that a Service Visit is overdue. • Contact Millipore SAS.

Continued on next page

Summary of Alert messages, Continued

Alert messages
(continued)
(continued)

TAP FEED CONDUCTIVITY > SP	<ul style="list-style-type: none"> • The Tap Water conductivity is below the set point. • Contact Millipore SAS.
THE NETWORK CABLE IS UNPLUGGED	<ul style="list-style-type: none"> • Check the Ethernet Cable plugged into the System and the computer. • Restart the System.
THIS IP ADDRESS IS ALREADY USED BY ANOTHER SYSTEM	<ul style="list-style-type: none"> • Contact your network administrator. • Restart the System.
TOC < 1 PPB	<ul style="list-style-type: none"> • Invalid TOC measurement. • Perform a TOC Curve Check. • Contact Millipore SAS.
TOC FEEDWATER < 15.3 MΩ.cm	<ul style="list-style-type: none"> • Due to low resistivity at the inlet of the UV Lamp, the TOC measurement can not be performed properly anymore. • Please replace the Q-Pak® Pack to ensure valid TOC indications.

Ordering Information

Consumables, Accessories and Systems

Milli-Q® Direct System

Item	Catalogue Number
Milli-Q® Direct 8	ZR0Q00800
Milli-Q® Direct 16	ZR0Q01600

Consumables

Item	Catalogue Number
Progard® T3 Cartridge	PROG000T3
Q-Pak® TEX Pack	QPAK00TEX
Q-Pak® TIX Pack	QPAK00TIX
BioPak® Ultrafilter	CDUFBI001
Millipak Express® 40 Final Filter	MPGP04001
EDS-Pak® Final Filter	EDSPAK001
VOC-Pak™ Final Filter	VOCPAK001
EDS-Pak® Installation Kit - ordered 1 time only for multiple EDS-Pak® uses.	EDSKIT001
UV 185 nm Lamp	ZMQUVLP01
UV 254 nm Lamp	ZLXUVLP01
ASM (Automatic Sanitisation Module) UV Lamp	ZLXUVLPL1

Accessories

Item	Catalogue Number
Reservoir 30 Litre	TANKPE030
Reservoir 60 Litre	TANKPE060
Reservoir 100 Litre	TANKPE100
Remote POD	ZMQSP0D02
ASM (Automatic Sanitisation Module) for Reservoir	TANKASMIN
Cabinet Wall Mounting Bracket	WMBSMT002
Footswitch (for Remote POD)	ZMQSFTS01
Remote POD Wall Mounting Bracket	WMBQPOD01
Water Sensor	ZFWATDET4

Note

Regularly scheduled preventive maintenance/calibration will help you obtain the best performance from your Millipore SAS water purification system throughout its entire lifetime.

Please contact your Millipore SAS representative to find the best options for your system including our maintenance programs.



FTPF11486



V3.0